



Biomet 3i Certain Internal	Dentsply Xive/Friadent	Nobel Biocare Replace (Tri-Lobe)
Biomet 3i Osseotite External	MegaGen AnyRidge	Osstem / Hiossen
Astra EV	MIS Seven	Straumann Bone Level
Astra OsseoSpeed	Multiunit Abutment	Straumann BLX
BioHorizons Internal	Neodent GM	Straumann Tissue Level
Camlog	Nobel Biocare Active	Zimmer TSV (Tapered Screw-Vent)
Conelog	Nobel Biocare Branemark	

Category	System	Kit	Reference Info
DESS AngleBase® 0.045 US			Engaging AngleBase: 37.XXX Non Engaging AngleBase: 36.XXX ** E: Engaging / NE: Non Engaging ** GH: Gingival Height

	Biomet 3i Certain Internal Desktop Scan Body: 50.014 , 50.015 , 50.016 IOS Scan Body: 52.014 , 52.015 , 52.016		
	DESS 3i® CER AngleBase® 0.045 v1.4		
	3i®_CER_WP_E 37.016	5.0 (6/5) Yellow	
	3i®_CER_RP_E 37.015	4.1 (4.0, 5/4) Blue	

	Biomet 3i Osseotite External Desktop Scan Body: 50.011 , 50.012 , 50.013 IOS Scan Body: 52.011 , 52.012 , 52.013		
	DESS 3i® OSS AngleBase® 0.045 v1.4		
	3i®_OSS_RP_E 37.012/3	4.1 (3.75, 4.0, 5/4) GH 3.0 mm	
	3i®_OSS_NP_E 37.011/3	3.4 (3.25, 4/3) GH 3.0 mm	
	3i®_OSS_RP_NE 36.012/3	4.1 (3.75, 4.0, 5/4) GH 3.0 mm	
	3i®_OSS_NP_NE 36.011/3	3.4 (3.25, 4/3) GH 3.0 mm	
	3i®_OSS_RP_E 37.012/2	4.1 (3.75, 4.0, 5/4) GH 1.5 mm	
	3i®_OSS_NP_E 37.011/2	3.4 (3.25, 4/3) GH 1.5 mm	
	3i®_OSS_RP_NE 36.012/2	4.1 (3.75, 4.0, 5/4) GH 1.5 mm	
	3i®_OSS_NP_NE 36.011/2	3.4 (3.25, 4/3) GH 1.5 mm	
	3i®_OSS_WP_E 37.013	5.0 (6/5) GH 0.3 mm	
	3i®_OSS_RP_E 37.012	4.1 (3.75, 4.0, 5/4) GH 0.3 mm	
	3i®_OSS_NP_E 37.011	3.4 (3.25, 4/3) GH 0.3 mm	
	3i®_OSS_WP_NE 36.013	5.0 (6/5) GH 0.3 mm	

	Astra EV Desktop Scan Body: 50.059 , 50.060 , 50.061 , 50.058 IOS Scan Body: 52.406 , 52.059 , 52.060 , 52.061 , 52.407 , 52.058		
	DESS AST EVO AngleBase® 0.045 v2.1		
	AST_EVO_4.8_NE_36.061	EV 4.8 Blue	
	AST_EVO_4.2_NE_36.060	EV 4.2 Yellow	
	AST_EVO_3.6_NE_36.059	EV 3.6 Purple	
	AST_EVO_4.8_E_37.061	EV 4.8 Blue	
	AST_EVO_4.2_E_37.060	EV 4.2 Yellow	
	AST_EVO_3.6_E_37.059	EV 3.6 Purple	
	AST_EVO_UNIABUT_36.058	On EV UniAbutment	

	Astra OsseoSpeed Desktop Scan Body: 50.024 , 50.025 IOS Scan Body: 52.023 , 52.024 , 52.025		
	DESS AST OSS AngleBase® 0.045 v1.3		
	AST_OSS_Lilac_E_37.025	4.5/5.0 Lilac	
	AST_OSS_Aqua_E_37.024	3.5/4.0 Aqua	
	AST_OSS_Lilac_NE_36.025	4.5/5.0 Lilac	

DESS AngleBase 0.045 US



Biomet 3i Certain Internal	Dentsply Xive/Friadent	Nobel Biocare Replace (Tri-Lobe)
Biomet 3i Osseotite External	MegaGen AnyRidge	Osstem / Hiossen
Astra EV	MIS Seven	Straumann Bone Level
Astra OsseoSpeed	Multiunit Abutment	Straumann BLX
BioHorizons Internal	Neodent GM	Straumann Tissue Level
Camlog	Nobel Biocare Active	Zimmer TSV (Tapered Screw-Vent)
Conelog	Nobel Biocare Branemark	

Category	System	Kit	Reference Info
DESS AngleBase® 0.045 US			Engaging AngleBase: 37.XXX Non Engaging AngleBase: 36.XXX ** E: Engaging / NE: Non Engaging ** GH: Gingival Height

 Zirkon Category: DESS Angled Base 0.045 US System: DESS BH INT AngleBase® 0.045 Kit: BH INTERNAL NP NE 36.087 BH INTERNAL NP NE 36.087 BH INTERNAL RP NE 36.088 BH INTERNAL WP NE 36.089 BH INTERNAL NP E 37.087 BH INTERNAL RP E 37.088 BH INTERNAL WP E 37.089	BioHorizons Internal		
	Desktop Scan Body: N/A		
	IOS Scan Body: 52.099 , 52.087 , 52.088 , 52.089		
	DESS BH INT AngleBase® 0.045	BH INTERNAL NP NE 36.087	3.5 Yellow
	v1.3	BH INTERNAL RP NE 36.088	4.5 Green
		BH INTERNAL WP NE 36.089	5.7 Blue
		BH INTERNAL NP E 37.087	3.5 Yellow
	BH INTERNAL RP E 37.088	4.5 Green	
	BH INTERNAL WP E 37.089	5.7 Blue	

 Zirkon Category: DESS Angled Base 0.045 US System: DESS CAMLOG AngleBase® 0.045 Kit: CAMLOG 4.3 NE 36.066 CAMLOG 4.3 NE 36.066 CAMLOG 4.3 E 37.066 CAMLOG 3.8 NE 36.065 CAMLOG 3.8 E 37.065	Camlog		
	Desktop Scan Body: 50.065 , 50.066		
	IOS Scan Body: 52.064 , 52.065 , 52.066 , 52.067		
	DESS CAMLOG AngleBase® 0.045	CAMLOG 4.3 NE 36.066	4.3 Pink
	v1.5	CAMLOG 4.3 E 37.066	4.3 Pink
		CAMLOG 3.8 NE 36.065	3.8 Yellow
	CAMLOG 3.8 E 37.065	3.8 Yellow	

 Zirkon Category: DESS Angled Base 0.045 US System: DESS CONELOG AngleBase® 0.0 Kit: CONELOG NP NE 36.074 CONELOG NP NE 36.074 CONELOG RP NE 36.075 CONELOG WP NE 36.076 CONELOG NP E 37.074 CONELOG RP E 37.075 CONELOG WP E 37.076	Conelog		
	Desktop Scan Body: N/A		
	IOS Scan Body: 52.074 , 52.075 , 52.076		
	DESS CONELOG AngleBase® 0.045	CONELOG NP NE 36.074	3.3 Silver
	v1.5	CONELOG RP NE 36.075	3.8 Yellow / 4.3 Pink
		CONELOG WP NE 36.076	5.0 Blue
		CONELOG NP E 37.074	3.3 Silver
	CONELOG RP E 37.075	3.8 Yellow / 4.3 Pink	
	CONELOG WP E 37.076	5.0 Blue	

 Zirkon Category: DESS Angled Base 0.045 US System: DESS DEN XIV AngleBase® 0.045 Kit: DEN_XIV_WP_E 37.040 DEN_XIV_WP_E 37.040 DEN_XIV_RP_E 37.039 DEN_XIV_NP_E 37.038 DEN_XIV_WP_NE 36.040 DEN_XIV_RP_NE 36.039 DEN_XIV_NP_NE 36.038	Dentsply Xive/Friadent		
	Desktop Scan Body: 50.038 , 50.039 , 50.040		
	IOS Scan Body: 52.038 , 52.039 , 52.040		
	DESS DEN XIV AngleBase® 0.045	DEN_XIV_WP_E 37.040	4.5 Blue
	v1.3	DEN_XIV_RP_E 37.039	3.8 Yellow
		DEN_XIV_NP_E 37.038	3.4 Silver
		DEN_XIV_WP_NE 36.040	4.5 Blue
		DEN_XIV_RP_NE 36.039	3.8 Yellow
	DEN_XIV_NP_NE 36.038	3.4 Silver	

 Zirkon Category: DESS Angled Base 0.045 US System: DESS MGN ANY AngleBase® 0.045 Kit: MGN ANY NE 36.057 MGN ANY NE 36.057 MGN ANY E 37.057	MegaGen AnyRidge		
	Desktop Scan Body: 50.057		
	IOS Scan Body: 52.057		
DESS MGN ANY AngleBase® 0.045	MGN ANY NE 36.057		
v1.3	MGN ANY E 37.057		

DESS AngleBase 0.045 US



Biomet 3i Certain Internal	Dentsply Xive/Friadent	Nobel Biocare Replace (Tri-Lobe)
Biomet 3i Osseotite External	MegaGen AnyRidge	Osstem / Hiossen
Astra EV	MIS Seven	Straumann Bone Level
Astra OsseoSpeed	Multiunit Abutment	Straumann BLX
BioHorizons Internal	Neodent GM	Straumann Tissue Level
Camlog	Nobel Biocare Active	Zimmer TSV (Tapered Screw-Vent)
Conelog	Nobel Biocare Branemark	

Category	System	Kit	Reference Info
DESS AngleBase® 0.045 US			Engaging AngleBase: 37.XXX Non Engaging AngleBase: 36.XXX ** E: Engaging / NE: Non Engaging ** GH: Gingival Height

	MIS Seven		
	Desktop Scan Body: 50.056 , 50.017 , 50.018		
	IOS Scan Body: 52.056 , 52.017 , 52.018		
	DESS MIS SEV AngleBase® 0.045	MIS_SEV_WP_NE_36.018	WP (5.0 & 6.0) Green
	v1.4	MIS_SEV_WP_E_37.018	WP (5.0 & 6.0) Green

	Multiunit Abutment		
	Desktop Scan Body: 50.007 , 50.047		
	IOS Scan Body: 52.007 , 52.047		
DESS MULTI-UNIT AngleBase® 0.045	MULTI-UNIT®_WP_NE_36.054	only for Nobel Branemark WP MUA	
v1.3	MULTI-UNIT®_RP_NE_36.007	DESS MUA RP components	

	Neodent GM (Grand Morse)		
	Desktop Scan Body: N/A		
	IOS Scan Body: 52.035		
DESS NEO GM AngleBase® 0.045	NEO_GM_NE_36.035_3	GH 2.5 mm	
v1.4	NEO_GM_E_37.035_3	GH 2.5 mm	

	Nobel Biocare Active / Replace Conical Connection (CC)		
	Desktop Scan Body: 50.068 , 50.041 , 50.042 , 50.069		
	IOS Scan Body: 52.068 , 52.041 , 52.042 , 52.069		
	DESS NOB ACT AngleBase® 0.045	NOB_ACT_RP_E_37.042	RP 4.3/5.0 Yellow
	v1.4	NOB_ACT_NP_E_37.041	NP 3.5 Pink

	Nobel Biocare Branemark (External)		
	Desktop Scan Body: 50.001 , 50.002 , 50.003		
	IOS Scan Body: 52.001 , 52.002 , 52.003		
	DESS NOB BRK AngleBase® 0.045	NOB_BRK_RP_E_37.002/3	RP 3.75/4.0 GH 2.5 mm
	v1.4	NOB_BRK_NP_E_37.001/3	NP 3.3 GH 2.5 mm
		NOB_BRK_RP_NE_36.002/3	RP 3.75/4.0 GH 2.5 mm
		NOB_BRK_NP_NE_36.001/3	NP 3.3 GH 2.5 mm
		NOB_BRK_RP_E_37.002/2	RP 3.75/4.0 GH 1.5 mm
		NOB_BRK_NP_E_37.001/2	NP 3.3 GH 1.5 mm
		NOB_BRK_RP_NE_36.002/2	RP 3.75/4.0 GH 1.5 mm
		NOB_BRK_NP_NE_36.001/2	NP 3.3 GH 1.5 mm
		NOB_BRK_RP_E_37.002	RP 3.75/4.0 GH 0.3 mm
		NOB_BRK_NP_E_37.001	NP 3.3 GH 0.3 mm
		NOB_BRK_RP_NE_36.002	RP 3.75/4.0 GH 0.3 mm
	NOB_BRK_NP_NE_36.001	NP 3.3 GH 0.3 mm	



Biomet 3i Certain Internal	Dentsply Xive/Friadent	Nobel Biocare Replace (Tri-Lobe)
Biomet 3i Osseotite External	MegaGen AnyRidge	Osstem / Hiossen
Astra EV	MIS Seven	Straumann Bone Level
Astra OsseoSpeed	Multiunit Abutment	Straumann BLX
BioHorizons Internal	Neodent GM	Straumann Tissue Level
Camlog	Nobel Biocare Active	Zimmer TSV (Tapered Screw-Vent)
Conelog	Nobel Biocare Branemark	

Category	System	Kit	Reference Info
DESS AngleBase® 0.045 US			Engaging AngleBase: 37.XXX Non Engaging AngleBase: 36.XXX ** E: Engaging / NE: Non Engaging ** GH: Gingival Height

<p>Zirkon</p> <p>Category: DESS Angled Base 0.045 US</p> <p>System: DESS NOB RPL AngleBase® 0.045</p> <p>Kit: NOB_RPL_WP_E_37.006</p> <p>NOB_RPL_WP_E_37.006 NOB_RPL_RP_E_37.005 NOB_RPL_NP_E_37.004 NOB_RPL_WP_NE_36.006 NOB_RPL_RP_NE_36.005 NOB_RPL_NP_NE_36.004</p>	Nobel Biocare Replace (Tri-Lobe) Desktop Scan Body: 50.004 , 50.005 , 50.006 , 50.097 IOS Scan Body: 52.004 , 52.005 , 52.006 , 52.097		
	DESS NOB RPL AngleBase® 0.045		
	v1.3	NOB_RPL_WP_E_37.006	WP 5.0 Blue
		NOB_RPL_RP_NE_36.005	RP 4.3 Yellow
		NOB_RPL_NP_E_37.004	NP 3.5 Pink
		NOB_RPL_WP_NE_36.006	WP 5.0 Blue
		NOB_RPL_RP_E_37.005	RP 4.3 Yellow
	NOB_RPL_NP_NE_36.004	NP 3.5 Pink	

<p>Zirkon</p> <p>Category: DESS Angled Base 0.045 US</p> <p>System: DESS OSS TS AngleBase® 0.045</p> <p>Kit: OSS_TS_Mini_NE_36.072</p> <p>OSS_TS_Mini_NE_36.072 OSS_TS_Mini_E_37.072 OSS_TS_STD_NE_36.073 OSS_TS_STD_E_37.073</p>	Osstem / Hiossen Desktop Scan Body: 50.072 , 50.073 IOS Scan Body: 52.072 , 52.073		
	DESS OSS TS AngleBase® 0.045		
	v2.0	OSS_TS_Mini_NE_36.072	Mini
		OSS_TS_Mini_E_37.072	Mini
		OSS_TS_STD_NE_36.073	Regular
		OSS_TS_STD_E_37.073	Regular

<p>Zirkon</p> <p>Category: DESS Angled Base 0.045 US</p> <p>System: DESS STR BLV AngleBase® 0.045</p> <p>Kit: STR_BLV_RC_E_37.044</p> <p>STR_BLV_RC_E_37.044 STR_BLV_NC_E_37.043 STR_BLV_RC_NE_36.044 STR_BLV_NC_NE_36.043</p>	Straumann Bone Level Desktop Scan Body: 50.043 , 50.044 IOS Scan Body: 52.043 , 52.044		
	DESS STR BLV AngleBase® 0.045		
	v1.4	STR_BLV_RC_E_37.044	RC 4.1/4.8 Pink
		STR_BLV_NC_E_37.043	NC 3.3 Yellow
		STR_BLV_RC_NE_36.044	RC 4.1/4.8 Pink
		STR_BLV_NC_NE_36.043	NC 3.3 Yellow

<p>Zirkon</p> <p>Category: DESS Angled Base 0.045 US</p> <p>System: DESS STR BLX AngleBase® 0.045</p> <p>Kit: STR_BLV_WB_NE_36.063</p> <p>STR_BLV_WB_NE_36.063 STR_BLV_WB_E_37.063 STR_BLV_RB_NE_36.062 STR_BLV_RB_E_37.062</p>	Straumann BLX Desktop Scan Body: N/A IOS Scan Body: 52.062		
	DESS STR BLX AngleBase® 0.045		
	v1.2	STR_BLV_WB_NE_36.063	5.0 Magenta, 5.5 Brown, 6.5 Black
		STR_BLV_WB_E_37.063	5.0 Magenta, 5.5 Brown, 6.5 Black
		STR_BLV_RB_NE_36.062	3.75 Red, 4.0 Grey, 4.5 Green
		STR_BLV_RB_E_37.062	3.75 Red, 4.0 Grey, 4.5 Green

<p>Zirkon</p> <p>Category: DESS Angled Base 0.045 US</p> <p>System: DESS STR STLX AngleBase® 0.045</p> <p>Kit: STR_STLV_WN_NE_36.047</p> <p>STR_STLV_WN_NE_36.047 STR_STLV_WN_E_37.045 STR_STLV_RN_E_37.009 STR_STLV_RN_NE_36.046</p>	Straumann Tissue Level Desktop Scan Body: 50.070 , 50.009 , 50.045 IOS Scan Body: 52.070 , 52.009 , 52.045		
	DESS STR STLX AngleBase® 0.045		
	v1.7	STR_STLV_WN_NE_36.047	WN 6.5
		STR_STLV_WN_E_37.045	WN 6.5
		STR_STLV_RN_E_37.009	RN 4.8
		STR_STLV_RN_NE_36.046	RN 4.8

DESS AngleBase 0.045 US



Biomet 3i Certain Internal	Dentsply Xive/Friadent	Nobel Biocare Replace (Tri-Lobe)
Biomet 3i Osseotite External	MegaGen AnyRidge	Osstem / Hiossen
Astra EV	MIS Seven	Straumann Bone Level
Astra OsseoSpeed	Multiunit Abutment	Straumann BLX
BioHorizons Internal	Neodent GM	Straumann Tissue Level
Camlog	Nobel Biocare Active	Zimmer TSV (Tapered Screw-Vent)
Conelog	Nobel Biocare Branemark	

Category	System	Kit	Reference Info												
DESS AngleBase® 0.045 US			Engaging AngleBase: 37.XXX Non Engaging AngleBase: 36.XXX ** E: Engaging / NE: Non Engaging ** GH: Gingival Height												
<p>Zirkon</p> Category: DESS Angled Base 0.045 US System: DESS ZIM SCV AngleBase® 0.04: Kit: ZIM_SCV_NP_E_37.017 ZIM_SCV_NP_E_37.017 ZIM_SCV_NP_NE_36.017 ZIM_SCV_RP_E_37.018 ZIM_SCV_RP_NE_36.018 ZIM_SCV_WP_E_37.019 ZIM_SCV_WP_NE_36.019	Zimmer TSV (Tapered Screw-Vent) Desktop Scan Body: 50.017 , 50.018 , 50.019 IOS Scan Body: 52.017 , 52.018 , 52.019														
	DESS ZIM SCV AngleBase® 0.045 v1.3		<table border="1"> <tr> <td>ZIM_SCV_NP E_37.017</td> <td>3.5 Green</td> </tr> <tr> <td>ZIM_SCV_NP NE_36.017</td> <td>3.5 Green</td> </tr> <tr> <td>ZIM_SCV_RP E_37.018</td> <td>4.5 Purple</td> </tr> <tr> <td>ZIM_SCV_RP NE_36.018</td> <td>4.5 Purple</td> </tr> <tr> <td>ZIM_SCV_WP E_37.019</td> <td>5.7 Yellow</td> </tr> <tr> <td>ZIM_SCV_WP_NE_36.019</td> <td>5.7 Yellow</td> </tr> </table>	ZIM_SCV_NP E_37.017	3.5 Green	ZIM_SCV_NP NE_36.017	3.5 Green	ZIM_SCV_RP E_37.018	4.5 Purple	ZIM_SCV_RP NE_36.018	4.5 Purple	ZIM_SCV_WP E_37.019	5.7 Yellow	ZIM_SCV_WP_NE_36.019	5.7 Yellow
	ZIM_SCV_NP E_37.017	3.5 Green													
	ZIM_SCV_NP NE_36.017	3.5 Green													
	ZIM_SCV_RP E_37.018	4.5 Purple													
	ZIM_SCV_RP NE_36.018	4.5 Purple													
	ZIM_SCV_WP E_37.019	5.7 Yellow													
	ZIM_SCV_WP_NE_36.019	5.7 Yellow													