

DESS Direct (2)



Astra EV	Straumann Bone Level on SRA (Screw-retained Abutment)
Astra OsseoSpeed	Straumann Tissue Level on Straumann Converter Abutment
DESS Multiunit Abutment with screw 19.006, 19.031, 19.306, & 19.331	DESS Multiunit Abutment with screw 19.018 & 19.406-P10

Category	System	Kit	Reference Info
DESS Direct			

<p>DESS UNIABUTMENT DIRECT (Uniabutment Scan Abutment)</p>	<p>Astra EV</p> <p>Desktop Scan Body: 50.058 IOS Scan Body: 52.058</p>		
	<p>DESS AST EVO Uniabutment DIRECT (Uniabutment scan) v1.6</p>	EV UNIABUTMENT	On EV UniAbutment
	<p>Astra OsseoSpeed</p> <p>Desktop Scan Body: 50.036, 50.037 IOS Scan Body: 52.036, 52.037</p>		
	<p>DESS AST Uniabutment® DIRECT (Uniabutment® Scan Abutment) v5.6</p>	20 45	On 20° UniAbutment On 45° UniAbutment

<p>DESS Multi-Unit DIRECT (Multi-unit Scan Abutment)</p>	<p>DESS Multiunit Abutment with screw 19.006, 19.031, 19.306, & 19.331</p> <p>Desktop Scan Body: 50.007, 50.047 IOS Scan Body: 52.007, 52.047</p>		
	<p>DESS Multi-Unit DIRECT (Multi-unit Scan Abutment) v5.5</p>	Multi-Unit	MUA Component info
		Multi-Unit WP	only for Nobel Branemark WP MUA

<p>DESS STR BLV SRA DIRECT (SRA Scan Abutment)</p>	<p>Straumann Bone Level on SRA (Screw-retained Abutment)</p> <p>Desktop Scan Body: N/A IOS Scan Body: 52.098</p>		
	<p>DESS STR BLV SRA DIRECT (SRA Scan Abutment) v1.2</p>	SRA Direct	On Straumann SRA (Screw-retained Abutment)

<p>DESS STR SYN DIRECT (Bars and bridges)</p>	<p>Straumann Tissue Level on Straumann Converter Abutment</p> <p>Desktop Scan Body: 50.010, 50.048 IOS Scan Body: 52.010, 52.048</p>		
	<p>DESS STR SYN DIRECT (Bars and bridges) v5.6</p>	STR SYN WN_51.048	WN 6.5 on Straumann Converter Abutment
		STR SYN RN_51.010	RN 4.8 on Straumann Converter Abutment

<p>DESS Multi-Unit Temporary</p>	<p>DESS Multiunit Abutment with screw 19.018 & 19.406-P10</p> <p>Desktop Scan Body: 50.007, 50.047 IOS Scan Body: 52.007, 52.047</p>		
	<p>DESS Multi-Unit Temporary v1.1</p>	Multi-Unit	MUA Component info

** screw for DESS Multi-Unit DIRECT Multi-Unit: **19.006** (UG), **19.055** (HEX 0.048"/1.20 mm), **19.065** (HEX 0.048"/1.22 mm), **19.110** (HEX 0.050"/1.27 mm) / Angled: **19.306** (TORX® Ball)
Multi-Unit WP: **19.031** (UG) / Angled: 19.331 (TORX® Ball)

** screw for DESS Multi-Unit Temporary Multi-Unit: **19.018** (HEX 0.050" / 1.27 mm) / Angled: **19.406-P101** (TORX® Ball)