

Abutment

Ti

Category

+

System

Kit

Category	System	Kit	Reference Info
DESS PRE-MILLED ABUTMENT CoCr			Premilled Blank Size: 71.XXX - 10 mmD 72.XXX - 14 mmD ** E: Engaging / NE: Non Engaging
<b>Biomet 3i Certain Internal</b>			
Desktop Scan Body: <a href="#">50.014</a> , <a href="#">50.015</a> , <a href="#">50.016</a> IOS Scan Body: <a href="#">52.014</a> , <a href="#">52.015</a> , <a href="#">52.016</a>			
	<b>DESS 3i® CER PRE-MILLED CoCr</b> v2.7	3i_CER_NP_E_71.014	3.4 (3.25, 4/3) Purple
		3i_CER_RP_E_71.015	4.1 (4.0, 5/4) Blue
		3i_CER_WP_E_71.016	5.0 (6/5) Yellow
<b>Biomet 3i Osseotite External</b>			
Desktop Scan Body: <a href="#">50.011</a> , <a href="#">50.012</a> , <a href="#">50.013</a> IOS Scan Body: <a href="#">52.011</a> , <a href="#">52.012</a> , <a href="#">52.013</a>			
	<b>DESS 3i® OSS PRE-MILLED CoCr</b> v2.2	3i_OSS_WP E 71.013	5.0 (6/5)
		3i_OSS_NP E 71.011	3.4 (3.25, 4/3)
		3i_OSS_RP E 71.012	4.1 (3.75, 4.0, 5/4)
<b>Astra EV</b>			
Desktop Scan Body: <a href="#">50.059</a> , <a href="#">50.060</a> , <a href="#">50.061</a> IOS Scan Body: <a href="#">52.406</a> , <a href="#">52.059</a> , <a href="#">52.060</a> , <a href="#">52.061</a> , <a href="#">52.407</a>			
	<b>DESS AST EVO PRE-MILLED CoCr</b> v2.8	AST EVO 3.6 E 71.059	EV 3.6 Purple
		AST EVO 4.2 E 71.060	EV 4.2 Yellow
		AST EVO 4.8 E 71.061	EV 4.8 Blue
		AST EVO 5.4 E 71.407	EV 5.4 Brown
<b>Astra Osseespeed</b>			
Desktop Scan Body: <a href="#">50.024</a> , <a href="#">50.025</a> IOS Scan Body: <a href="#">52.023</a> , <a href="#">52.024</a> , <a href="#">52.025</a>			
	<b>DESS AST OSS PRE-MILLED CoCr</b> v1.3	AST OSS Aqua E 71.024	3.5/4.0 Aqua
		AST OSS Lilac E 71.025	4.5/5.0 Lilac
<b>BioHorizons External</b>			
Desktop Scan Body: <a href="#">50.020</a> , <a href="#">50.021</a> , <a href="#">50.022</a> IOS Scan Body: <a href="#">52.020</a> , <a href="#">52.021</a> , <a href="#">52.022</a>			
	<b>DESS BH EXTERNAL PRE-MILLED CoCr</b> v2.7	BH EXTERNAL RP E 71.021	4.0 (4.2) Green
		BH EXTERNAL NP E 71.020	3.5 (3.7) Yellow
		BH EXTERNAL WP E 71.022	5.0 (5.2) Blue
<b>Camlog</b>			
Desktop Scan Body: <a href="#">50.065</a> , <a href="#">50.066</a> IOS Scan Body: <a href="#">52.064</a> , <a href="#">52.065</a> , <a href="#">52.066</a> , <a href="#">52.067</a>			
	<b>DESS CAMLOG PRE-MILLED CoCr</b> v2.4	CAMLOG 3.8 E 71.065	3.8 Yellow
		CAMLOG 4.3 E 71.066	4.3 Pink
		CAMLOG 3.3 E 71.064	3.3 Silver
		CAMLOG 5.0 E 71.067	5.0 Blue
<b>Conelog</b>			
Desktop Scan Body: N/A IOS Scan Body: <a href="#">52.074</a> , <a href="#">52.075</a> , <a href="#">52.076</a>			
	<b>DESS CONELOG PRE-MILLED CoCr</b> v2.3	CONELOG NP 3.3 71.074	3.3 Silver
		CONELOG RP 3.8 71.075	3.8 Yellow / 4.3 Pink
		CONELOG WP 5.0 71.076	5.0 Blue

Abutment	Ti		+
	Category	<input type="text"/>	
	System	<input type="text"/>	
	Kit	<input type="text"/>	

Category	System	Kit	Reference Info
DESS PRE-MILLED ABUTMENT CoCr			Premilled Blank Size: 71.XXX - 10 mmD 72.XXX - 14 mmD ** E: Engaging / NE: Non Engaging
<b>Dentsply Ankylos</b>			
Desktop Scan Body: <a href="#">50.049</a> IOS Scan Body: <a href="#">52.049</a>			
DESS DEN ANK PRE-MILLED CoCr v2.3		DEN_ANK_DIRECT 71.250	Ankylos/X
<b>Dentsply Xive/Friadent</b>			
Desktop Scan Body: <a href="#">50.038</a> , <a href="#">50.039</a> , <a href="#">50.040</a> IOS Scan Body: <a href="#">52.038</a> , <a href="#">52.039</a> , <a href="#">52.040</a>			
DESS DEN XIV PRE-MILLED CoCr v2.3		DEN_XIV_NP_E 71.038	3.4 Silver
		DEN_XIV_RP_E 71.039	3.8 Yellow
		DEN_XIV_WP_E 71.040	4.5 Blue
		DEN_XIV_5.5_E 71.071	5.5 Red
<b>MegaGen AnyRidge</b>			
Desktop Scan Body: <a href="#">50.057</a> IOS Scan Body: <a href="#">52.057</a>			
DESS MGEN ANY PRE-MILLED CoCr v2.3		MGEN ANY E 71.057	
<b>MIS Seven</b>			
Desktop Scan Body: <a href="#">50.056</a> , <a href="#">50.017</a> , <a href="#">50.018</a> IOS Scan Body: <a href="#">52.056</a> , <a href="#">52.017</a> , <a href="#">52.018</a>			
DESS MIS SEV PRE-MILLED CoCr v1.5		MIS SEV SP E 71.095	SP (3.75 & 4.2) Purple
		MIS SEV WP E 71.096	WP (5.0 & 6.0) Green
<b>Neodent GM (Grand Morse)</b>			
Desktop Scan Body: N/A IOS Scan Body: <a href="#">52.035</a>			
DESS NEO GM PRE-MILLED CoCr v1.1		NEO GM E 71.035	
<b>Nobel Biocare Active / Replace Conical Connection (CC)</b>			
Desktop Scan Body: <a href="#">50.068</a> , <a href="#">50.041</a> , <a href="#">50.042</a> , <a href="#">50.069</a> IOS Scan Body: <a href="#">52.068</a> , <a href="#">52.041</a> , <a href="#">52.042</a> , <a href="#">52.069</a>			
DESS NOB ACT PRE-MILLED CoCr v2.7		NOB ACT RP E 71.042	RP 4.3/5.0 Yellow
		NOB ACT NP E 71.041	NP 3.5 Pink
		NOB ACT 3.0 E 71.068	3.0 Silver
		NOB ACT WP E 71.069	WP 5.5 Blue
<b>Nobel Biocare Branemark (External)</b>			
Desktop Scan Body: <a href="#">50.001</a> , <a href="#">50.002</a> , <a href="#">50.003</a> IOS Scan Body: <a href="#">52.001</a> , <a href="#">52.002</a> , <a href="#">52.003</a>			
DESS NOB BRK PRE-MILLED CoCr v2.3		NOB BRK RP E 71.002	RP 3.75/4.0
		NOB BRK WP E 71.003	WP 5.0
		NOB BRK NP E 71.001	NP 3.3

Abutment	Ti	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Category	<input type="text"/>			
	System	<input type="text"/>			
	Kit	<input type="text"/>			

Category	System	Kit	Reference Info
DESS PRE-MILLED ABUTMENT CoCr			Premilled Blank Size: 71.XXX - 10 mmD 72.XXX - 14 mmD ** E: Engaging / NE: Non Engaging
	<b>Nobel Biocare Replace (Tri-Lobe)</b>		
	Desktop Scan Body: <a href="#">50.004</a> , <a href="#">50.005</a> , <a href="#">50.006</a> , <a href="#">50.097</a> IOS Scan Body: <a href="#">52.004</a> , <a href="#">52.005</a> , <a href="#">52.006</a> , <a href="#">52.097</a>		
<b>DESS NOB RPL PRE-MILLED CoCr</b> v2.5	NOB RPL NP E 71.004	NP 3.5 Pink	
	NOB RPL RP E 71.005	RP 4.3 Yellow	
	NOB RPL WP E 71.006	WP 5.0 Blue	
	NOB RPL 6.0 E 71.097	6.0 Green	
	<b>Osstem / Hiossen</b>		
	Desktop Scan Body: <a href="#">50.072</a> , <a href="#">50.073</a> IOS Scan Body: <a href="#">52.072</a> , <a href="#">52.073</a>		
<b>DESS OSS TS PRE-MILLED CoCr</b> v2.5	OSS TS Mini NP E 71.072	Mini	
	OSS TS Std RP E 71.073	Regular	
	<b>Straumann Bone Level</b>		
	Desktop Scan Body: <a href="#">50.043</a> , <a href="#">50.044</a> IOS Scan Body: <a href="#">52.043</a> , <a href="#">52.044</a>		
<b>DESS STR BLV PRE-MILLED CoCr</b> v2.5	STR BLV NC E 71.043	NC 3.3 Yellow	
	STR BLV RC E 71.044	RC 4.1/4.8 Pink	
	<b>Straumann BLX</b>		
	Desktop Scan Body: N/A IOS Scan Body: <a href="#">52.062</a>		
<b>DESS STR BLX PRE-MILLED CoCr</b> v1.1	STR BLX WB E 71.063	5.0 Magenta, 5.5 Brown, 6.5 Black	
	<b>Straumann Tissue Level</b>		
	Desktop Scan Body: <a href="#">50.070</a> , <a href="#">50.009</a> , <a href="#">50.045</a> IOS Scan Body: <a href="#">52.070</a> , <a href="#">52.009</a> , <a href="#">52.045</a>		
<b>DESS STR STLV PRE-MILLED CoCr</b> v2.8	STR STLV WN 71.045	WN 6.5	
	STR STLV RN 71.009	RN 4.8	
	STR STLV NNC 71.070	NNC 3.5	
	<b>Zimmer TSV (Tapered Screw-Vent)</b>		
	Desktop Scan Body: <a href="#">50.017</a> , <a href="#">50.018</a> , <a href="#">50.019</a> IOS Scan Body: <a href="#">52.017</a> , <a href="#">52.018</a> , <a href="#">52.019</a>		
<b>DESS ZIM SCV PRE-MILLED CoCr</b> v2.5	ZIM SCV NP E 71.017	3.5 Green	
	ZIM SCV RP E 71.018	4.5 Purple	
	ZIM SCV WP E 71.019	5.7 Yellow	