

DESS Ti Base 0.045



Biomet 3i Certain Internal	Conelog	Nobel Biocare Active
Biomet 3i Osseotite External	Dentsply Ankylos	Nobel Biocare Branemark
Astra EV	Dentsply Xive/Friadent	Nobel Biocare Replace (Tri-Lobe)
Astra OsseoSpeed	MegaGen AnyRidge	Osstem / Hiossen
BioHorizons External	MIS Seven	Straumann Bone Level
BioHorizons Internal	Multiunit Abutment	Straumann BLX
Camlog	Neodent GM	Straumann Tissue Level
		Straumann Tissue Level on Straumann Converter Abutment
		Zimmer TSV (Tapered Screw-Vent)

Category	System	Kit	Reference Info
DESS Ti Base 0.045			Engaging Ti Base: 16.XXX Non Engaging Ti Base: 15.XXX ** E: Engaging / NE: Non Engaging ** GH: Gingival Height

 Zirkon Category: DESS Ti Base 0.045 System: DESS 3i® CER Ti Base 0.045 Kit: 3i CER NP E 16.214	Biomet 3i Certain Internal Desktop Scan Body: 50.014, 50.015, 50.016 IOS Scan Body: 52.014, 52.015, 52.016		
	DESS 3i® CER Ti Base 0.045 v8.4	3i CER NP E 16.214 3i CER RP E 16.215 3i CER RP E 16.315 3i CER WP E 16.216 3i CER WP E 16.316 3i CER NP NE 15.214 3i CER NP NE 15.314 3i CER RP NE 15.215 3i CER RP NE 15.315 3i CER WP NE 15.216 3i CER WP NE 15.316 3i CER NP E 16.014 3i CER NP NE 15.014 3i CER WP E 16.016 3i CER WP NE 15.016 3i CER RP E 16.015 3i CER RP NE 15.015 3i CER NP E 16.314	3.4 (3.25, 4/3) Purple GH 1.5 mm 4.1 (4.0, 5/4) Blue GH 1.5 mm 4.1 (4.0, 5/4) Blue GH 3.0 mm 5.0 (6/5) Yellow GH 1.5 mm 5.0 (6/5) Yellow GH 3.0 mm 3.4 (3.25, 4/3) Purple GH 1.5 mm 3.4 (3.25, 4/3) Purple GH 3.0 mm 4.1 (4.0, 5/4) Blue GH 1.5 mm 4.1 (4.0, 5/4) Blue GH 3.0 mm 5.0 (6/5) Yellow GH 1.5 mm 5.0 (6/5) Yellow GH 3.0 mm 3.4 (3.25, 4/3) Purple GH 0.3 mm 3.4 (3.25, 4/3) Purple GH 0.3 mm 5.0 (6/5) Yellow GH 0.3 mm 5.0 (6/5) Yellow GH 0.3 mm 4.1 (4.0, 5/4) Blue GH 0.3 mm 4.1 (4.0, 5/4) Blue GH 0.3 mm 3.4 (3.25, 4/3) Purple GH 3.0 mm

 Zirkon Category: DESS Ti Base 0.045 System: DESS 3i® OSS Ti Base 0.045 Kit: 3i OSS NP E 16.011	Biomet 3i Osseotite External Desktop Scan Body: 50.011, 50.012, 50.013 IOS Scan Body: 52.011, 52.012, 52.013		
	DESS 3i® OSS Ti Base 0.045 v8.5	3i OSS NP E 16.011 3i OSS NP NE 15.011 3i OSS RP E 16.012 3i OSS RP NE 15.012 3i OSS WP E 16.013 3i OSS WP NE 15.013 3i OSS NP NE 15.211 3i OSS NP NE 15.311 3i OSS RP NE 15.212 3i OSS RP NE 15.312 3i OSS WP NE 15.213 3i OSS WP NE 15.313 3i OSS NP E 16.211 3i OSS NP E 16.311 3i OSS RP E 16.212 3i OSS RP E 16.312 3i OSS WP E 16.213 3i OSS WP E 16.313	3.4 (3.25, 4/3) GH 0.3 mm 3.4 (3.25, 4/3) GH 0.3 mm 4.1 (3.75, 4.0, 5/4) GH 0.3 mm 4.1 (3.75, 4.0, 5/4) GH 0.3 mm 5.0 (6/5) GH 0.3 mm 5.0 (6/5) GH 0.3 mm 3.4 (3.25, 4/3) GH 1.5 mm 3.4 (3.25, 4/3) GH 3.0 mm 4.1 (3.75, 4.0, 5/4) GH 1.5 mm 4.1 (3.75, 4.0, 5/4) GH 3.0 mm 5.0 (6/5) GH 1.5 mm 5.0 (6/5) GH 3.0 mm 3.4 (3.25, 4/3) GH 1.5 mm 3.4 (3.25, 4/3) GH 3.0 mm 4.1 (3.75, 4.0, 5/4) GH 1.5 mm 4.1 (3.75, 4.0, 5/4) GH 3.0 mm 5.0 (6/5) GH 1.5 mm 5.0 (6/5) GH 3.0 mm

 Zirkon Category: DESS Ti Base 0.045 System: DESS AST EVO Ti Base 0.045 Kit: AST EVO 3.6 E 16.059	Astra EV Desktop Scan Body: 50.059, 50.060, 50.061, 50.058 IOS Scan Body: 52.406, 52.059, 52.060, 52.061, 52.407, 52.058		
	DESS AST EVO Ti Base 0.045 v4.3	AST EVO 3.6 E 16.059 AST EVO 4.2 E 16.060 AST EVO 4.8 E 16.061 AST EVO 3.6 NE 15.359 AST EVO 4.2 NE 15.360 AST EVO 4.8 NE 15.361 AST EVO 3.6 E 16.359 AST EVO 4.2 E 16.360 AST EVO 4.8 E 16.361 AST EVO UNIABUTMENT NE 15.058 AST EVO 3.6 NE 15.059 AST EVO 4.2 NE 15.060 AST EVO 4.8 NE 15.061 AST EVO 5.4 E 16.407 AST EVO 5.4 NE 15.407	EV 3.6 Purple GH 1.0 mm EV 4.2 Yellow GH 1.0 mm EV 4.8 Blue GH 1.0 mm EV 3.6 Purple GH 3.0 mm EV 4.2 Yellow GH 3.0 mm EV 4.8 Blue GH 3.0 mm EV 3.6 Purple GH 3.0 mm EV 4.2 Yellow GH 3.0 mm EV 4.8 Blue GH 3.0 mm On EV UniAbutment EV 3.6 Purple GH 1.0 mm EV 4.2 Yellow GH 1.0 mm EV 4.8 Blue GH 1.0 mm EV 5.4 Brown GH 3.0 mm EV 5.4 Brown GH 3.0 mm

DESS Ti Base 0.045



Biomet 3i Certain Internal	Conelog	Nobel Biocare Active
Biomet 3i Osseotite External	Dentsply Ankylos	Nobel Biocare Branemark
Astra EV	Dentsply Xive/Friadent	Nobel Biocare Replace (Tri-Lobe)
Astra OsseoSpeed	MegaGen AnyRidge	Osstem / Hiossen
BioHorizons External	MIS Seven	Straumann Bone Level
BioHorizons Internal	Multiunit Abutment	Straumann BLX
Camlog	Neodent GM	Straumann Tissue Level
		Straumann Tissue Level on Straumann Converter Abutment
		Zimmer TSV (Tapered Screw-Vent)

Category	System	Kit	Reference Info
DESS Ti Base 0.045			Engaging Ti Base: 16.XXX Non Engaging Ti Base: 15.XXX ** E: Engaging / NE: Non Engaging ** GH: Gingival Height

Zirkon Category: DESS Ti Base 0.045 System: DESS AST OSS Ti Base 0.045 Kit: AST OSS Yellow NE 15.023 AST OSS Yellow NE 15.023 AST_OSS_Yellow_E 16.023 AST_OSS_Lilac_E 16.325 AST_OSS_Aqua E 16.024 AST_OSS_Aqua NE 15.024 AST_OSS_Lilac_E 16.025 AST_OSS_UNIABUTMENT_20 NE 15.036 AST_OSS_UNIABUTMENT_45 NE 15.037 AST_OSS_Lilac_NE 15.025 AST_OSS_Lilac NE 15.325 AST_OSS_Aqua E 16.324 AST_OSS_Aqua NE 15.324 AST_OSS_Yellow_NE 15.023/2 AST_OSS_Yellow_NE 15.023/3 AST_OSS_Yellow_NE 16.023/2 AST_OSS_Yellow_NE 16.023/3 AST_OSS_Aqua E 16.024/2 AST_OSS_Aqua E 15.024/2 AST_OSS_Lilac_E 16.025/2 AST_OSS_Lilac_E 15.025/2	Astra OsseoSpeed Desktop Scan Body: 50.024, 50.025, 50.036, 50.037 IOS Scan Body: 52.023, 52.024, 52.025, 52.036, 52.037		
	DESS AST OSS Ti Base 0.045 v8.9	AST_OSS_Yellow_NE 15.023 3.0 Yellow GH 1.0 mm AST_OSS_Yellow_E 16.023 3.0 Yellow GH 1.0 mm AST_OSS_Lilac_E 16.325 4.5/5.0 Lilac GH 3.0 mm AST_OSS_Aqua E 16.024 3.5/4.0 Aqua GH 1.1 mm AST_OSS_Aqua NE 15.024 3.5/4.0 Aqua GH 1.1 mm AST_OSS_Lilac_E 16.025 4.5/5.0 Lilac GH 0.9 mm AST_OSS_UNIABUTMENT_20 NE 15.036 On 20° UniAbutment AST_OSS_UNIABUTMENT_45 NE 15.037 On 45° UniAbutment AST_OSS_Lilac_NE 15.025 4.5/5.0 Lilac GH 0.9 mm AST_OSS_Lilac NE 15.325 4.5/5.0 Lilac GH 3.0 mm AST_OSS_Aqua E 16.324 3.5/4.0 Aqua GH 3.0 mm AST_OSS_Aqua NE 15.324 3.5/4.0 Aqua GH 3.0 mm AST_OSS_Yellow_NE 15.023/2 3.0 Yellow GH 2.0 mm AST_OSS_Yellow_NE 15.023/3 3.0 Yellow GH 3.0 mm AST_OSS_Yellow_NE 16.023/2 3.0 Yellow GH 2.0 mm AST_OSS_Yellow_NE 16.023/3 3.0 Yellow GH 3.0 mm AST_OSS_Aqua E 16.024/2 3.5/4.0 Aqua GH 2.0 mm AST_OSS_Aqua E 15.024/2 3.5/4.0 Aqua GH 2.0 mm AST_OSS_Lilac_E 16.025/2 4.5/5.0 Lilac GH 2.0 mm AST_OSS_Lilac_E 15.025/2 4.5/5.0 Lilac GH 2.0 mm	

Zirkon Category: DESS Ti Base 0.045 System: DESS BH EXT Ti Base 0.045 Kit: BH_EXTERNAL NP E 16.051 BH_EXTERNAL NP E 16.051 BH_EXTERNAL_NP NE 15.051 BH_EXTERNAL_RP E 16.052 BH_EXTERNAL_RP NE 15.052 BH_EXTERNAL_WP E 16.053 BH_EXTERNAL_WP NE 15.053	BioHorizons External Desktop Scan Body: 50.020, 50.021, 50.022 IOS Scan Body: 52.020, 52.021, 52.022		
	DESS BH EXT Ti Base 0.045 v6.7	BH_EXTERNAL NP E 16.051 3.5 Yellow GH 0.3 mm BH_EXTERNAL_NP NE 15.051 3.5 Yellow GH 0.3 mm BH_EXTERNAL_RP E 16.052 4.0 Green GH 0.3 mm BH_EXTERNAL_RP NE 15.052 4.0 Green GH 0.3 mm BH_EXTERNAL_WP E 16.053 5.0 Blue GH 0.3 mm BH_EXTERNAL_WP NE 15.053 5.0 Blue GH 0.3 mm	

Zirkon Category: DESS Ti Base 0.045 System: DESS BH INT Ti Base 0.045 Kit: BH_INTERNAL NP NE 15.087 BH_INTERNAL NP NE 15.087 BH_INTERNAL_RP NE 15.088 BH_INTERNAL_WP NE 15.089 BH_INTERNAL NP E 16.087 BH_INTERNAL_RP E 16.088 BH_INTERNAL_WP E 16.089	BioHorizons Internal Desktop Scan Body: N/A IOS Scan Body: 52.099, 52.087, 52.088, 52.089		
	DESS BH INT Ti Base 0.045 v41.4	BH_INTERNAL NP NE 15.087 3.5 Yellow GH 0.8 mm BH_INTERNAL_RP NE 15.088 4.5 Green GH 0.3 mm BH_INTERNAL_WP NE 15.089 5.7 Blue GH 0.3 mm BH_INTERNAL NP E 16.087 3.5 Yellow GH 0.8 mm BH_INTERNAL_RP E 16.088 4.5 Green GH 0.3 mm BH_INTERNAL_WP E 16.089 5.7 Blue GH 0.3 mm	

Zirkon Category: DESS Ti Base 0.045 System: DESS CAMLOG Ti Base 0.045 Kit: CAM_3.8 E 16.065 CAM_3.8 E 16.065 CAM_3.8_E 16.265 CAM_3.8_E 16.365 CAM_3.8 NE 15.065 CAM_3.8_NE 15.265 CAM_3.8_NE 15.365 CAM_4.3 E 16.066 CAM_4.3_E 16.266 CAM_4.3_NE 15.066 CAM_4.3_NE 15.266 CAM_4.3_NE 15.366 CAM_3.3_NE 15.064 CAM_5.0_NE 15.067 CAM_3.3_E 16.064 CAM_5.0_E 16.067	Camlog Desktop Scan Body: 50.065, 50.066 IOS Scan Body: 52.064, 52.065, 52.066, 52.067		
	DESS CAMLOG Ti Base 0.045 v3.8	CAM_3.8 E 16.065 3.8 Yellow GH 0.5 mm CAM_3.8_E 16.265 3.8 Yellow GH 1.5 mm CAM_3.8_E 16.365 3.8 Yellow GH 3.0 mm CAM_3.8 NE 15.065 3.8 Yellow GH 0.5 mm CAM_3.8_NE 15.265 3.8 Yellow GH 1.5 mm CAM_3.8_NE 15.365 3.8 Yellow GH 3.0 mm CAM_4.3 E 16.066 4.3 Pink GH 0.3 mm CAM_4.3_E 16.266 4.3 Pink GH 1.5 mm CAM_4.3_E 16.366 4.3 Pink GH 3.0 mm CAM_4.3_NE 15.066 4.3 Pink GH 0.3 mm CAM_4.3_NE 15.266 4.3 Pink GH 1.5 mm CAM_4.3_NE 15.366 4.3 Pink GH 3.0 mm CAM_3.3_NE 15.064 3.3 Sliver GH 0.5 mm CAM_5.0_NE 15.067 5.0 Blue GH 0.3 mm CAM_3.3_E 16.064 3.3 Sliver GH 0.5 mm CAM_5.0_E 16.067 5.0 Blue GH 0.3 mm	

DESS Ti Base 0.045



Biomet 3i Certain Internal	Conelog	Nobel Biocare Active
Biomet 3i Osseotite External	Dentsply Ankylos	Nobel Biocare Branemark
Astra EV	Dentsply Xive/Friadent	Nobel Biocare Replace (Tri-Lobe)
Astra OsseoSpeed	MegaGen AnyRidge	Osstem / Hiossen
BioHorizons External	MIS Seven	Straumann Bone Level
BioHorizons Internal	Multiunit Abutment	Straumann BLX
Camlog	Neodent GM	Straumann Tissue Level
		Straumann Tissue Level on Straumann Converter Abutment
		Zimmer TSV (Tapered Screw-Vent)

Category	System	Kit	Reference Info
DESS Ti Base 0.045			Engaging Ti Base: 16.XXX Non Engaging Ti Base: 15.XXX ** E: Engaging / NE: Non Engaging ** GH: Gingival Height

	Conelog Desktop Scan Body: N/A IOS Scan Body: 52.074, 52.075, 52.076		
	DESS CONELOG Ti Base 0.045 v1.5	CONELOG NP NE 15.074 CONELOG RP NE 15.075 CONELOG WP NE 15.076 CONELOG NP E 16.074 CONELOG RP E 16.075 CONELOG WP E 16.076 CONELOG NP NE 15.074/3 CONELOG RP NE 15.075/3 CONELOG WP NE 15.076/3 CONELOG NP E 16.074/3 CONELOG RP E 16.075/3 CONELOG WP E 16.076/3	3.3 Silver GH 1.4 mm 3.8 Yellow / 4.3 Pink GH 1.4 mm 5.0 Blue GH 1.4 mm 3.3 Silver GH 1.4 mm 3.8 Yellow / 4.3 Pink GH 1.4 mm 5.0 Blue GH 1.4 mm 3.3 Silver GH 3.0 mm 3.8 Yellow / 4.3 Pink GH 3.0 mm 5.0 Blue GH 3.0 mm 3.3 Silver GH 3.0 mm 3.8 Yellow / 4.3 Pink GH 3.0 mm 5.0 Blue GH 3.0 mm

	Dentsply Ankylos Desktop Scan Body: 50.049 IOS Scan Body: 52.049		
	DESS DEN ANK Ti Base 0.045 v6.5	DEN ANK NE GH1.5 15.050 DEN ANK E GH1.5 16.050 DEN ANK E GH2.5 16.250 DEN ANK NE GH2.5 15.250	Ankylos/C GH 1.5 mm Ankylos/X GH 1.5 mm Ankylos/X GH 2.5 mm Ankylos/C GH 2.5 mm

	Dentsply Xive/Friadent Desktop Scan Body: 50.038, 50.039, 50.040 IOS Scan Body: 52.038, 52.039, 52.040		
	DESS DEN XIV Ti Base 0.045 v6.5	DEN XIV 5.5 NE 15.071 DEN XIV 3.4 NE 15.238 DEN XIV 3.4 NE 15.338 DEN XIV 3.8 NE 15.239 DEN XIV 3.8 NE 15.339 DEN XIV 4.5 NE 15.240 DEN XIV 4.5 NE 15.340 DEN XIV 3.4 E 16.238 DEN XIV 3.4 E 16.338 DEN XIV 3.8 E 16.239 DEN XIV 3.8 E 16.339 DEN XIV 4.5 E 16.240 DEN XIV 4.5 E 16.340 DEN XIV 5.5 E 16.071 DEN XIV 3.4 E 16.038 DEN XIV 3.4 NE 15.038 DEN XIV 3.8 E 16.039 DEN XIV 3.8 NE 15.039 DEN XIV 4.5 E 16.040 DEN XIV 4.5 NE 15.040 DEN XIV 5.5 NE 15.271 DEN XIV 5.5 NE 15.371 DEN XIV 5.5 E 16.271 DEN XIV 5.5 E 16.371	5.5 Red GH 0.3 mm 3.4 Silver GH 1.5 mm 3.4 Silver GH 3.0 mm 3.8 Yellow GH 1.5 mm 3.8 Yellow GH 3.0 mm 4.5 Blue GH 1.5 mm 4.5 Blue GH 3.0 mm 3.4 Silver GH 1.5 mm 3.4 Silver GH 3.0 mm 3.8 Yellow GH 1.5 mm 3.8 Yellow GH 3.0 mm 4.5 Blue GH 1.5 mm 4.5 Blue GH 3.0 mm 5.5 Red GH 0.3 mm 3.4 Silver GH 0.3 mm 3.4 Silver GH 0.3 mm 3.8 Yellow GH 0.3 mm 3.8 Yellow GH 0.3 mm 4.5 Blue GH 0.3 mm 4.5 Blue GH 0.3 mm 5.5 Red GH 1.5 mm 5.5 Red GH 3.0 mm 5.5 Red GH 1.5 mm 5.5 Red GH 3.0 mm

	MegaGen AnyRidge Desktop Scan Body: 50.057 IOS Scan Body: 52.057		
	DESS MGN ANY Ti Base 0.045 v2.7	MGN ANY E 16.057 MGN ANY E 16.357 MGN ANY NE 15.057 MGN ANY NE 15.357	GH 2.0 mm GH 3.0 mm GH 2.0 mm GH 3.0 mm

DESS Ti Base 0.045



Biomet 3i Certain Internal	Conelog	Nobel Biocare Active
Biomet 3i Osseotite External	Dentsply Ankylos	Nobel Biocare Branemark
Astra EV	Dentsply Xive/Friadent	Nobel Biocare Replace (Tri-Lobe)
Astra OsseoSpeed	MegaGen AnyRidge	Osstem / Hiossen
BioHorizons External	MIS Seven	Straumann Bone Level
BioHorizons Internal	Multiunit Abutment	Straumann BLX
Camlog	Neodent GM	Straumann Tissue Level
		Zimmer TSV (Tapered Screw-Vent)

Category	System	Kit	Reference Info
DESS Ti Base 0.045			Engaging Ti Base: 16.XXX Non Engaging Ti Base: 15.XXX ** E: Engaging / NE: Non Engaging ** GH: Gingival Height

 Zirkon Category: DESS Ti Base 0.045 System: DESS MIS SEV Ti Base 0.045 Kit: MIS_SEV_SP E 16.095 MIS_SEV_SP NE 15.095 MIS_SEV_WP E 16.096 MIS_SEV_WP NE 15.295 MIS_SEV_SP NE 15.395 MIS_SEV_SP E 16.295 MIS_SEV_SP E 16.395 MIS_SEV_WP E 16.296 MIS_SEV_WP E 16.396 MIS_SEV_WP NE 15.296 MIS_SEV_WP NE 15.396 MIS_SEV_3.0 E 16.056 MIS_SEV_3.0 NE 15.056 MIS_SEV_WP NE 15.096 MIS_SEV_3.0 E 16.256 MIS_SEV_3.0 E 16.356 MIS_SEV_3.0 NE 15.256 MIS_SEV_3.0 NE 15.356	MIS Seven Desktop Scan Body: 50.056, 50.017, 50.018 IOS Scan Body: 52.056, 52.017, 52.018		
	DESS MIS SEV Ti Base 0.045 v6.8	MIS_SEV_SP E 16.095 SP (3.75 & 4.2) Purple GH 0.8 mm MIS_SEV_SP NE 15.095 SP (3.75 & 4.2) Purple GH 0.8 mm MIS_SEV_WP E 16.096 WP (5.0 & 6.0) Green GH 0.3 mm MIS_SEV_SP NE 15.295 SP (3.75 & 4.2) Purple GH 1.5 mm MIS_SEV_SP NE 15.395 SP (3.75 & 4.2) Purple GH 3.0 mm MIS_SEV_SP E 16.295 SP (3.75 & 4.2) Purple GH 1.5 mm MIS_SEV_SP E 16.395 SP (3.75 & 4.2) Purple GH 3.0 mm MIS_SEV_WP E 16.296 WP (5.0 & 6.0) Green GH 1.5 mm MIS_SEV_WP E 16.396 WP (5.0 & 6.0) Green GH 3.0 mm MIS_SEV_WP NE 15.296 WP (5.0 & 6.0) Green GH 1.5 mm MIS_SEV_WP NE 15.396 WP (5.0 & 6.0) Green GH 3.0 mm MIS_SEV_3.0 E 16.056 NP (3.3) Yellow GH 0.8 mm MIS_SEV_3.0 NE 15.056 NP (3.3) Yellow GH 0.8 mm MIS_SEV_WP NE 15.096 WP (5.0 & 6.0) Green GH 0.3 mm MIS_SEV_3.0 E 16.256 NP (3.3) Yellow GH 1.5 mm MIS_SEV_3.0 E 16.356 NP (3.3) Yellow GH 3.0 mm MIS_SEV_3.0 NE 15.256 NP (3.3) Yellow GH 1.5 mm MIS_SEV_3.0 NE 15.356 NP (3.3) Yellow GH 3.0 mm	

 Zirkon Category: DESS Ti Base 0.045 System: DESS MULTI-UNIT Ti Base 0.045 Kit: Multi_Unit_RP_NE 15.007/8 Multi_Unit_RP NE 15.007/8 Multi_Unit_RP NE 15.007 Multi_Unit_WP NE 15.054	Multiunit Abutment Desktop Scan Body: 50.007, 50.047 IOS Scan Body: 52.007, 52.047		
	DESS MULTI-UNIT Ti Base 0.045 v6.7	Multi_Unit_RP_NE 15.007/8 Chimney Height 8 mm DESS MUA Multi_Unit_RP NE 15.007 Chimney Height 6 mm DESS MUA Multi_Unit_WP NE 15.054 only for Nobel Branemark WP MUA	

 Zirkon Category: DESS Ti Base 0.045 System: DESS NEO GM Ti Base 0.045 Kit: NEO_GM_NE 15.035 NEO_GM_NE 16.035 NEO_GM_NE 15.235 NEO_GM_E 16.235 NEO_GM_NE 15.335 NEO_GM_E 16.335	Neodent GM (Grand Morse) Desktop Scan Body: N/A IOS Scan Body: 52.035		
	DESS NEO GM Ti Base 0.045 v1.3	NEO_GM_NE 15.035 GH 0.8 mm NEO_GM_E 16.035 GH 0.8 mm NEO_GM_NE 15.235 GH 1.5 mm NEO_GM_E 16.235 GH 1.5 mm NEO_GM_NE 15.335 GH 2.5 mm NEO_GM_E 16.335 GH 2.5 mm	

 Zirkon Category: DESS Ti Base 0.045 System: DESS NOB ACT Ti Base 0.045 Kit: NOB_ACT_3.0 E 16.068 NOB_ACT_3.0 E 16.068 NOB_ACT_3.0 NE 15.068 NOB_ACT_WP NE 15.069 NOB_ACT_NP NE 15.341 NOB_ACT_RP NE 15.342 NOB_ACT_NP NE 16.341 NOB_ACT_RP E 16.342 NOB_ACT_NP NE 15.041 NOB_ACT_RP E 16.042 NOB_ACT_RP NE 15.042 NOB_ACT_WP E 16.069 NOB_ACT_NP E 16.041 NOB_ACT_NP E 16.041/2 NOB_ACT_RP E 16.042/2 NOB_ACT_NP NE 15.041/2 NOB_ACT_RP NE 15.042/2	Nobel Biocare Active / Replace Conical Connection (CC) Desktop Scan Body: 50.068, 50.041, 50.042, 50.069 IOS Scan Body: 52.068, 52.041, 52.042, 52.069		
	DESS NOB ACT Ti Base 0.045 v6.2	NOB_ACT_3.0 E 16.068 3.0 Silver GH 1.5 mm NOB_ACT_3.0 NE 15.068 3.0 Silver GH 1.5 mm NOB_ACT_WP NE 15.069 WP 5.5 Blue GH 1.5 mm NOB_ACT_NP NE 15.341 NP 3.5 Pink GH 3.0 mm NOB_ACT_RP NE 15.342 RP 4.3/5.0 Yellow GH 3.0 mm NOB_ACT_NP E 16.341 NP 3.5 Pink GH 3.0 mm NOB_ACT_RP E 16.342 RP 4.3/5.0 Yellow GH 3.0 mm NOB_ACT_NP NE 15.041 NP 3.5 Pink GH 1.2 mm NOB_ACT_RP E 16.042 RP 4.3/5.0 Yellow GH 1.2 mm NOB_ACT_RP NE 15.042 RP 4.3/5.0 Yellow GH 1.2 mm NOB_ACT_WP E 16.069 WP 5.5 Blue GH 1.5 mm NOB_ACT_NP E 16.041 NP 3.5 Pink GH 1.2 mm NOB_ACT_NP E 16.041/2 NP 3.5 Pink GH 2.0 mm NOB_ACT_RP E 16.042/2 RP 4.3/5.0 Yellow GH 2.0 mm NOB_ACT_NP NE 15.041/2 NP 3.5 Pink GH 2.0 mm NOB_ACT_RP NE 15.042/2 RP 4.3/5.0 Yellow GH 2.0 mm	

DESS Ti Base 0.045



Biomet 3i Certain Internal	Conelog	Nobel Biocare Active
Biomet 3i Osseotite External	Dentsply Ankylos	Nobel Biocare Branemark
Astra EV	Dentsply Xive/Friadent	Nobel Biocare Replace (Tri-Lobe)
Astra OsseoSpeed	MegaGen AnyRidge	Osstem / Hiossen
BioHorizons External	MIS Seven	Straumann Bone Level
BioHorizons Internal	Multiunit Abutment	Straumann BLX
Camlog	Neodent GM	Straumann Tissue Level
		Straumann Tissue Level on Straumann Converter Abutment
		Zimmer TSV (Tapered Screw-Vent)

Category	System	Kit	Reference Info
DESS Ti Base 0.045			Engaging Ti Base: 16.XXX Non Engaging Ti Base: 15.XXX ** E: Engaging / NE: Non Engaging ** GH: Gingival Height

 Zirkon Category: DESS Ti Base 0.045 System: DESS NOB BRK Ti Base 0.045 Kit: NOB BRK NP E 16.201	Nobel Biocare Branemark (External) Desktop Scan Body: 50.001, 50.002, 50.003 IOS Scan Body: 52.001, 52.002, 52.003		
	DESS NOB BRK Ti Base 0.045 v6.5	NOB BRK NP E 16.201 NP 3.3 GH 1.5 mm NOB BRK NP E 16.301 NP 3.3 GH 3.0 mm NOB BRK RP E 16.202 RP 3.75/4.0 GH 1.5 mm NOB BRK RP E 16.302 RP 3.75/4.0 GH 3.0 mm NOB BRK WP E 16.203 WP 5.0 GH 1.5 mm NOB BRK WP E 16.303 WP 5.0 GH 3.0 mm NOB BRK NP NE 15.201 NP 3.3 GH 1.5 mm NOB BRK NP NE 15.301 NP 3.3 GH 3.0 mm NOB BRK RP NE 15.202 RP 3.75/4.0 GH 1.5 mm NOB BRK RP NE 15.302 RP 3.75/4.0 GH 3.0 mm NOB BRK WP NE 15.203 WP 5.0 GH 1.5 mm NOB BRK WP NE 15.303 WP 5.0 GH 3.0 mm NOB BRK NP E 16.001 NP 3.3 GH 0.3 mm NOB BRK NP NE 15.001 NP 3.3 GH 0.3 mm NOB BRK RP E 16.002 RP 3.75/4.0 GH 0.3 mm NOB BRK RP NE 15.002 RP 3.75/4.0 GH 0.3 mm NOB BRK WP E 16.003 WP 5.0 GH 0.3 mm NOB BRK WP NE 15.003 WP 5.0 GH 0.3 mm	

 Zirkon Category: DESS Ti Base 0.045 System: DESS NOB RPL Ti Base 0.045 Kit: NOB RPL NP NE 15.204	Nobel Biocare Replace (Tri-Lobe) Desktop Scan Body: 50.004, 50.005, 50.006, 50.097 IOS Scan Body: 52.004, 52.005, 52.006, 52.097		
	DESS NOB RPL Ti Base 0.045 v6.5	NOB RPL NP NE 15.204 NP 3.5 Pink GH 1.5 mm NOB RPL NP NE 15.304 NP 3.5 Pink GH 3.0 mm NOB RPL RP NE 15.205 RP 4.3 Yellow GH 1.5 mm NOB RPL RP NE 15.305 RP 4.3 Yellow GH 3.0 mm NOB RPL WP NE 15.206 WP 5.0 Blue GH 1.5 mm NOB RPL WP NE 15.306 WP 5.0 Blue GH 3.0 mm NOB RPL 6.0 NE 15.097 6.0 Green GH 0.3 mm NOB RPL NP E 16.204 NP 3.5 Pink GH 1.5 mm NOB RPL NP E 16.304 NP 3.5 Pink GH 3.0 mm NOB RPL RP E 16.205 RP 4.3 Yellow GH 1.5 mm NOB RPL RP E 16.305 RP 4.3 Yellow GH 3.0 mm NOB RPL WP E 16.206 WP 5.0 Blue GH 1.5 mm NOB RPL WP E 16.306 WP 5.0 Blue GH 3.0 mm NOB RPL 6.0 E 16.097 6.0 Green GH 0.3 mm NOB RPL NP E 16.004 NP 3.5 Pink GH 0.3 mm NOB RPL NP NE 15.004 NP 3.5 Pink GH 0.3 mm NOB RPL RP E 16.005 RP 4.3 Yellow GH 0.3 mm NOB RPL RP NE 15.005 RP 4.3 Yellow GH 0.3 mm NOB RPL WP E 16.006 WP 5.0 Blue GH 0.3 mm NOB RPL WP NE 15.006 WP 5.0 Blue GH 0.3 mm	

 Zirkon Category: DESS Ti Base 0.045 System: DESS OSS TS Ti Base 0.045 Kit: OSS TS Mini NE 15.072	Osstem / Hiossen Desktop Scan Body: 50.072, 50.073 IOS Scan Body: 52.072, 52.073		
	DESS OSS TS Ti Base 0.045 v2.5	OSS TS Mini NE 15.072 Mini GH 1.0 mm OSS TS Mini E 16.072 Mini GH 1.0 mm OSS TS STD NE 15.073 Regular GH 1.0 mm OSS TS STD E 16.073 Regular GH 1.0 mm OSS TS Mini NE 15.372 Mini GH 3.0 mm OSS TS Mini E 16.372 Mini GH 3.0 mm OSS TS STD NE 15.373 Regular GH 3.0 mm OSS TS STD E 16.373 Regular GH 3.0 mm	

DESS Ti Base 0.045



Biomet 3i Certain Internal	Conelog	Nobel Biocare Active
Biomet 3i Osseotite External	Dentsply Ankylos	Nobel Biocare Branemark
Astra EV	Dentsply Xive/Friadent	Nobel Biocare Replace (Tri-Lobe)
Astra OsseoSpeed	MegaGen AnyRidge	Osstem / Hiossen
BioHorizons External	MIS Seven	Straumann Bone Level
BioHorizons Internal	Multiunit Abutment	Straumann BLX
Camlog	Neodent GM	Straumann Tissue Level
		Straumann Tissue Level on Straumann Converter Abutment
		Zimmer TSV (Tapered Screw-Vent)

Category	System	Kit	Reference Info
DESS Ti Base 0.045			Engaging Ti Base: 16.XXX Non Engaging Ti Base: 15.XXX ** E: Engaging / NE: Non Engaging ** GH: Gingival Height

 Zirkon Category: DESS Ti Base 0.045 System: DESS STR BLV Ti Base 0.045 Kit: STR BLV on SRA NE 15.098 STR BLV on SRA NE 15.098 STR BLV on SRA E 16.098 STR_BLV_NC NE 15.343 STR_BLV_RC NE 15.344 STR_BLV_NC E 16.343 STR_BLV_NC E 16.043 STR_BLV_NC NE 15.043 STR_BLV_RC E 16.044 STR_BLV_RC NE 15.044 STR_BLV_RC E 16.344 STR_BLV_NC E 16.043/2 STR_BLV_NC NE 15.043/2 STR_BLV_RC E 16.044/2 STR_BLV_RC NE 15.044/2	Straumann Bone Level Desktop Scan Body: 50.043, 50.044 IOS Scan Body: 52.043, 52.044, 52.098		
	DESS STR BLV Ti Base 0.045 v7.1	STR BLV on SRA NE 15.098 STR BLV on SRA E 16.098 STR_BLV_NC NE 15.343 STR_BLV_RC NE 15.344 STR_BLV_NC E 16.343 STR_BLV_NC E 16.043 STR_BLV_NC NE 15.043 STR_BLV_RC E 16.044 STR_BLV_RC NE 15.044 STR_BLV_RC E 16.344 STR_BLV_NC E 16.043/2 STR_BLV_NC NE 15.043/2 STR_BLV_RC E 16.044/2 STR_BLV_RC NE 15.044/2	On Straumann SRA (Screw-retained Abutment) On Straumann SRA (Screw-retained Abutment) NC 3.3 Yellow GH 3.0 mm RC 4.1/4.8 Pink GH 3.0 mm NC 3.3 Yellow GH 3.0 mm NC 3.3 Yellow GH 1.2 mm NC 3.3 Yellow GH 1.2 mm RC 4.1/4.8 Pink GH 1.2 mm RC 4.1/4.8 Pink GH 1.2 mm RC 4.1/4.8 Pink GH 3.0 mm NC 3.3 Yellow GH 2.0 mm NC 3.3 Yellow GH 2.0 mm RC 4.1/4.8 Pink GH 2.0 mm RC 4.1/4.8 Pink GH 2.0 mm


 Zirkon Category: DESS Ti Base 0.045 System: DESS STR BLX Ti Base 0.045 Kit: STR BLX RB E 16.062/2 STR BLX RB E 16.062/2 STR BLX RB E 16.062/3 STR BLX RB E 16.062/4 STR BLX RB NE 15.062/2 STR BLX RB NE 15.062/3 STR BLX RB NE 15.062/4 STR BLX WB E 16.063 STR BLX WB E 16.063/2 STR BLX WB NE 15.063 STR BLX WB NE 15.063/2	Straumann BLX Desktop Scan Body: N/A IOS Scan Body: 52.062		
	DESS STR BLX Ti Base 0.045 v1.2	STR BLX RB E 16.062/2 STR BLX RB E 16.062/3 STR BLX RB E 16.062/4 STR BLX RB NE 15.062/2 STR BLX RB NE 15.062/3 STR BLX RB NE 15.062/4 STR BLX WB E 16.063 STR BLX WB E 16.063/2 STR BLX WB NE 15.063 STR BLX WB NE 15.063/2	3.75 Red, 4.0 Grey, 4.5 Green GH 1.5 mm 3.75 Red, 4.0 Grey, 4.5 Green GH 2.5 mm 3.75 Red, 4.0 Grey, 4.5 Green GH 3.5 mm 3.75 Red, 4.0 Grey, 4.5 Green GH 1.5 mm 3.75 Red, 4.0 Grey, 4.5 Green GH 2.5 mm 3.75 Red, 4.0 Grey, 4.5 Green GH 3.5 mm 5.0 Magenta, 5.5 Brown, 6.5 Black GH 0.8 mm 5.0 Magenta, 5.5 Brown, 6.5 Black GH 1.5 mm 5.0 Magenta, 5.5 Brown, 6.5 Black GH 0.8 mm 5.0 Magenta, 5.5 Brown, 6.5 Black GH 1.5 mm

 Zirkon Category: DESS Ti Base 0.045 System: DESS STR STL V Ti Base 0.045 Kit: STR STL V NNC NE 15.070 STR STL V NNC NE 15.070 STR STL V NNC E 16.070 STR STL V RN NE 15.046 STR STL V WN NE 15.045 STR STL V WN E 16.045 STR STL V RN E 16.009	Straumann Tissue Level Desktop Scan Body: 50.070, 50.009, 50.045 IOS Scan Body: 52.070, 52.009, 52.045		
	DESS STR STL V Ti Base 0.045 v7.0	STR STL V NNC NE 15.070 STR STL V NNC E 16.070 STR STL V RN NE 15.046 STR STL V WN NE 15.045 STR STL V WN E 16.045 STR STL V RN E 16.009	NNC 3.5 GH 0.5 mm NNC 3.5 GH 0.5 mm RN 4.8 GH 0.4 mm WN 6.5 GH 0.3 mm WN 6.5 GH 0.3 mm RN 4.8 GH 0.4 mm

 Zirkon Category: DESS Ti Base 0.045 System: DESS STR SYN Ti Base 0.045 Kit: STR SYN WN E 16.048 STR SYN WN E 16.048 STR SYN RN NE 15.010 STR SYN WN NE 15.048 STR SYN RN E 16.010	Straumann Tissue Level on Straumann Converter Abutment Desktop Scan Body: 50.010, 50.048 IOS Scan Body: 52.010, 52.048		
	DESS STR SYN Ti Base 0.045 v6.7	STR SYN WN E 16.048 STR SYN RN NE 15.010 STR SYN WN NE 15.048 STR SYN RN E 16.010	WN 6.5 GH 0.3 mm RN 4.8 GH 0.3 mm WN 6.5 GH 0.3 mm RN 4.8 GH 0.3 mm

DESS Ti Base 0.045



Biomet 3i Certain Internal	Conelog	Nobel Biocare Active	
Biomet 3i Osseotite External	Dentsply Ankylos	Nobel Biocare Branemark	
Astra EV	Dentsply Xive/Friadent	Nobel Biocare Replace (Tri-Lobe)	
Astra OsseoSpeed	MegaGen AnyRidge	Osstem / Hiossen	
BioHorizons External	MIS Seven	Straumann Bone Level	
BioHorizons Internal	Multiunit Abutment	Straumann BLX	Straumann Tissue Level on Straumann Converter Abutment
Camlog	Neodent GM	Straumann Tissue Level	Zimmer TSV (Tapered Screw-Vent)
Category	System	Kit	Reference Info
DESS Ti Base 0.045			Engaging Ti Base: 16.XXX Non Engaging Ti Base: 15.XXX ** E: Engaging / NE: Non Engaging ** GH: Gingival Height
 Zirkon Category: DESS Ti Base 0.045 System: DESS ZIM SCV Ti Base 0.045 Kit: ZIM_SCV_NP_NE_15.217 ZIM_SCV_NP_NE_15.217 ZIM_SCV_NP_NE_15.317 ZIM_SCV_RP_NE_15.218 ZIM_SCV_RP_NE_15.318 ZIM_SCV_WP_NE_15.219 ZIM_SCV_WP_NE_15.319 ZIM_SCV_NP_E_16.217 ZIM_SCV_NP_E_16.317 ZIM_SCV_RP_E_16.218 ZIM_SCV_RP_E_16.318 ZIM_SCV_WP_E_16.219 ZIM_SCV_WP_E_16.319 ZIM_SCV_NP_E_16.017 ZIM_SCV_NP_NE_15.017 ZIM_SCV_RP_E_16.018 ZIM_SCV_RP_NE_15.018 ZIM_SCV_WP_E_16.019 ZIM_SCV_WP_NE_15.019	Zimmer TSV (Tapered Screw-Vent) Desktop Scan Body: 50.017 , 50.018 , 50.019 IOS Scan Body: 52.017 , 52.018 , 52.019		
	DESS ZIM SCV Ti Base 0.045 v6.5	ZIM_SCV_NP_NE_15.217 3.5 Green GH 1.5 mm ZIM_SCV_NP_NE_15.317 3.5 Green GH 3.0 mm ZIM_SCV_RP_NE_15.218 4.5 Purple GH 1.5 mm ZIM_SCV_RP_NE_15.318 4.5 Purple GH 3.0 mm ZIM_SCV_WP_NE_15.219 5.7 Yellow GH 1.5 mm ZIM_SCV_WP_NE_15.319 5.7 Yellow GH 3.0 mm ZIM_SCV_NP_E_16.217 3.5 Green GH 1.5 mm ZIM_SCV_NP_E_16.317 3.5 Green GH 3.0 mm ZIM_SCV_RP_E_16.218 4.5 Purple GH 1.5 mm ZIM_SCV_RP_E_16.318 4.5 Purple GH 3.0 mm ZIM_SCV_WP_E_16.219 5.7 Yellow GH 1.5 mm ZIM_SCV_WP_E_16.319 5.7 Yellow GH 3.0 mm ZIM_SCV_NP_E_16.017 3.5 Green GH 0.8 mm ZIM_SCV_NP_NE_15.017 3.5 Green GH 0.8 mm ZIM_SCV_RP_E_16.018 4.5 Purple GH 0.3 mm ZIM_SCV_RP_NE_15.018 4.5 Purple GH 0.3 mm ZIM_SCV_WP_E_16.019 5.7 Yellow GH 0.3 mm ZIM_SCV_WP_NE_15.019 5.7 Yellow GH 0.3 mm	