



DESS

> DENTAL SMART SOLUTIONS



DESS

LIBRARY GUIDE

FOR HYBRID ABUTMENTS, BRIDGES, BARS AND
PROSTHESIS

Dental Wings library for implant and abutment level scan

For DWOS CAD system the library draws distinction between the following groups:



- DESS SYNOCTA® DIRECT (Bars and Bridges) - for multiple restorations connected directly with synOcta® abutment. The scanning needs to be performed with Scan Abutment for synOcta®: 41.50.010/41.52.010 (RN) or 41.50.048/41.52.048 (WN). **Note:** This connection is non- engaging only (See Table 1.2).
- DESS Ti-Base - for restorations connected with the implant via the Titanium Interface. The scanning needs to be performed with the Scan Body for the specific implant connection. **Note:** This category includes parts **scanned from abutment level** - Multi-Unit abutment, UniAbutment® and synOcta® abutment - with its specific Scan Abutment (see Table 1.1).
- DESS UNIABUTMENT® DIRECT (UniAbutment® Scan Abutment) - for restorations connected directly with UniAbutment®. The scanning needs to be performed with the Scan Abutment for UniAbutment® system – 41.50.036/41.52.036 (20°), 41.50.037/41.52.037 (45°) or 41.50.058/41.52.058. The desired UniAbutment® needs to be mounted physically on the implant analogue or model needs to be poured with UniAbutment® analogues (41.14.036, 41.14.037 or 41.14.058), (See Table 1.3).
- DESS Multi-Unit DIRECT (Multi-Unit Scan Abutment) - for restorations connected directly with Multi-Unit abutment. The scanning needs to be performed with the Scan Abutment for Multi-Unit system – 41.50.007/41.52.007 (RP) and 41.50.047/41.52.047 (WP). The desired Multi-Unit abutment needs to be mounted physically on the implant analogue or the model needs to be poured with Multi- unit® analogues (41.14.007 or 41.14.047), (See Table 1.2).
- DESS AURUMBase® - for restorations connected with the implant via the AURUMBase®. The scanning needs to be performed with the Scan Body for the specific implant connection. **Note:** This category includes parts **scanned from abutment level** - Multi- unit® abutment - with its specific Scan Abutment (see Table 1.3).
- DESS Pre-milled Ti – for restorations connected with the implant via the pre-milled abutment made in Titanium (see Table 1.6). The scanning needs to be performed with the Scan Body for the specific implant connection. **Note:** The maximum width of the piece is restricted to 10 mm.
- DESS Pre-milled CoCr – for restorations connected with the implant via the pre-milled abutment made in CoCr (see Table 1.7). The scanning needs to be performed with the Scan Body for the specific implant connection. **Note:** The maximum width of the piece is restricted to 14 mm.

Take care to select the correct DESS Scan Abutment and the corresponding library file, as otherwise the milled prosthesis can differ from desired outcome and errors can occur, for example the level of interface settlement.

Table 1.1 DESS TI-Base

	Platform	Analogue	Type	Kit	Scan Abutments	
Nobel Active®	3.0	14.068		NOB_ACT_3.0 NE 15.068	52.068	
					53.100	
				NOB_ACT_3.0 NE 16.068	52.068	
					52.200	
	NP/3.5	14.041		NOB_ACT_NP NE 15.041 15.341	50.041	52.041
					53.100	
				NOB_ACT_NP E 16.041 16.341	50.041	52.041
					52.200	
	RP/4.3	14.042		NOB_ACT_RP NE 15.042 15.342	50.042	52.042
					53.100	
				NOB_ACT_RP E 16.042 16.342	50.042	52.042
					52.200	
	WP/5.5	14.069		NOB_ACT_WP NE 15.069	52.069	
					53.100	
				NOB_ACT_WP NE 16.069	52.069	
					52.200	
Nobel Replace Select™	NP/3.5	14.004		NOB_RPL_NP NE 15.004 15.204 15.304	50.004	52.004
					53.100	
				NOB_RPL_NP E 16.004 16.204 16.304	50.004	52.004
					52.200	
	RP/4.3	14.005		NOB_RPL_RP NE 15.005 15.205 15.305	50.005	52.005
					53.100	
				NOB_RPL_RP E 16.005 16.205 16.305	50.005	52.005
					52.200	
	WP/5.0	14.006		NOB_RPL_WP NE 15.006 15.206 15.306	50.006	52.006
					53.100	
				NOB_RPL_WP E 16.006 16.206 16.306	50.006	52.006
					52.200	
	6.0	14.097		NOB_RPL_NP 6.0 15.097	52.097	
					53.100	
				NOB_RPL_NP 6.0 16.097	52.097	
					52.200	
Nobel Brånemark®	NP/3.5	14.001		NOB_BRK_NP NE 15.001 15.201 15.301	50.001	52.001
					53.100	
				NOB_BRK_NP E 16.001 16.201 16.301	50.001	52.001
					52.200	
	RP/4.1	14.002		NOB_BRK_RP NE 15.002 15.202 15.302	50.002	52.002
					53.100	
				NOB_BRK_RP E 16.002 16.202 16.302	50.002	52.002
					52.200	
	WP/5.0	14.003		NOB_BRK_WP NE 15.003 15.203 15.303	50.003	
					52.003	
				NOB_BRK_WP E 16.003 16.203 16.303	50.003	
					52.003	



non-engaging



non-engaging



Scan Abutment
for desktop scanning
50.0XX



Scan Abutment
for intra-oral scanning
52.0XX



Scan Abutment
on Ti base
53.100 53.200

	Platform	Analogue	Type	Kit	Scan Abutments
3i Osseotite®	NP/3.4	14.011		3I_OSS_NP NE 15.011	50.011 52.011
				15.211	53.100
				3I_OSS_NP E 16.011	50.011 52.011
	RP/4.1	14.012		15.212	52.200
				16.311	50.012 52.012
				3I_OSS_RP NE 15.312	53.100
	WP/5.0	14.013		16.011	50.012 52.012
				16.212	52.200
				3I_OSS_WP NE 15.013	50.013 52.013
				15.213	53.100
				16.313	50.013 52.013
				3I_OSS_WP E 16.213	52.200
3i Certain®	NP/3.4	14.014		3I_CER_NP NE 15.014	50.014 52.014
				15.214	53.100
				15.314	50.014 52.014
	RP/4.1	14.015		16.014	52.200
				16.214	50.015 52.015
				3I_CER_RP NE 15.015	53.100
	WP/5.0	14.016		15.315	50.015 52.015
				16.015	52.200
				3I_CER_RP E 16.215	50.016 52.016
				16.315	53.100
				15.016	50.016 52.016
				3I_CER_WP NE 15.216	52.200
Straumann® Soft Tissue Level	NNC/3.5	14.070		3I_CER_WP E 16.216	50.016 52.016
				16.316	52.200
				NOB_ACT_WP NE 15.069	52.069
	RN/4.8	14.009		NOB_ACT_WP NE 16.069	53.100
					52.069
				STR_STLV_RN NE 15.046	52.200
	RN on synOCTA®	14.010		STR_STLV_RN E 16.009	50.009
					52.009
				STR_SYN_RN NE 15.010	50.010
	WN/6.5	14.045		STR_SYN_RN E 16.010	52.010
					50.010
				STR_STLV_WP NE 15.045	50.045
	WN on synOCTA®	14.046		STR_STLV_WP E 16.045	52.045
					50.045
				STR_SYN_WN NE 15.048	50.015
				STR_SYN_WN E 16.048	52.015
					50.015
					52.015

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

 Scan Abutment
on Ti base
53.100 53.200

	Platform	Analogue	Type	Kit	Scan Abutments	
Straumann® Bone Level	NC/3.3	14.043		STR_BLV_NC NE 15.043 15.343	50.043	52.043 53.100
				TR_BLV_NC E 16.043 16.343	50.043	52.043 52.200
	RC/4.1	14.044		STR_BLV_RC NE 15.044 15.344	50.044	52.044 53.100
				TR_BLV_RC E 16.044 16.344	50.044	52.044 52.200
	SRA	14.098		STR_BLV_SRA 15.098	52.098	
				STR_BLV_SRA 16.098	52.098	
Zimmer Screw-Vent®	NP/3.5	14.017		ZIM_SCV_NP NE 15.017 15.317	50.017	52.017 53.100
				ZIM_SCV_NP E 16.017 16.217 16.317	50.017	52.017 52.200
	RP/4.5	14.018		ZIM_SCV_RP NE 15.018 15.318	50.018	52.018 53.100
				ZIM_SCV_RP E 16.018 16.218 16.318	50.018	52.018 52.200
	WP/5.7	14.019		ZIM_SCV_WP NE 15.019 15.219 15.319	50.019	52.019 53.100
				ZIM_SCV_WP E 16.019 16.219 16.319	50.019	52.019 52.200
Astra tech implant system™ EV	EV/3.6	14.059		AST_EVO 3.6_NE 15.059 15.359	50.059	52.059 53.100
				AST_EVO 3.6_E 16.059 16.359	50.059	52.059 52.200
	EV/4.2	14.060		AST_EVO 4.2_NE 15.060 15.360	50.060	52.060 53.100
				AST_EVO 4.2_E 16.060 16.360	50.060	52.060 52.200
	EV/4.8	14.061		AST_EVO 4.8_NE 15.061 15.361	50.061	52.061 53.100
				AST_EVO 4.8_E 16.061 16.361	50.061	52.061 52.200
	Uniabutment®	14.058		AST_EVO uniabutment_NE 15.058	50.058 52.058	
				AST_EVO uniabutment_NE 16.058	50.058 52.058	

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

 Scan Abutment
on Ti base
53.100 53.200

	Platform	Analogue	Type	Kit	Scan Abutments	
Astra tech osseospeed™	NP/3.0	14.023		AST_OSS_NP NE 15.023	52.023	
				AST_OSS_NP E 16.023	50.059	52.059
	RP/3.5-4.0	14.024		AST_OSS_RP NE 15.024 15.324	50.024	52.024
				AST_OSS_RP E 16.024 16.324	50.024	52.024
	WP/4.5-5.0	14.025		AST_OSS_WP NE 15.025 15.325	50.025	52.025
				AST_OSS_WP E 16.025 16.325	50.025	52.025
	Uniabutment® 20°	14.036		AST_OSS_Uniabutment® 20° NE 15.036	50.036	52.036
Uniabutment® 45°	14.037		AST_OSS_Uniabutment® 45° NE 15.037	50.037	52.037	
Friadent XIVE®	NP/3.4	14.038		DEN_XIV_NP NE 15.038	50.038	52.038
				DEN_XIV_NP E 16.038	50.038	52.038
	RP/3.8	14.039		DEN_XIV_RP NE 15.039	50.039	52.039
				DEN_XIV_RP E 16.039	50.039	52.039
	WP/4.5	14.040		DEN_XIV_WP NE 15.040	50.040	52.040
				DEN_XIV_WP E 16.040	50.040	52.040
	5.5	14.071		DEN_XIV_5.5 NE 15.071	52.071	
			DEN_XIV_5.4 E 16.071	52.071		
Ankylos® C/X	All platforms, height GH 1,5	14.048		DEN ANK NE GH1 15.050	50.049	52.049
		14.049		DEN ANK E GH1 16.050	50.049	52.049
	All platforms, height GH 2,5	14.048		DEN ANK NE GH2 15.250	50.049	52.049
		14.049		DEN ANK E GH2 16.250	50.049	52.049

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

 Scan Abutment
on Ti base
53.100 53.200

	Platform	Analogue	Type	Kit	Scan Abutments	
Osstem® TS	Mini Yellow	14.072		OSS_TS MINI_NE 15.072	50.072	52.072
				OSS_TS MINI_E 16.072	50.072	52.072
	STD Regular	14.073		OSS_TS STD_NE 15.372	53.100	52.200
				OSS_TS STD_E 16.372	50.073	52.073
Camlog®	3.3	14.064		CAM_3.3 NE 15.064	50.073	52.073
				CAM_3.3 E 16.064	53.100	52.200
	3.8	14.065		CAM_3.8 NE 15.065	50.065	52.065
				CAM_3.8 E 16.065	53.100	52.200
	4.3	14.066		CAM_4.3 NE 15.266	50.065	52.065
				CAM_4.3 E 16.266	50.066	52.066
	5.0	14.067		CAM_5.0 NE 15.067	53.100	52.200
				CAM_5.03 E 16.067	52.067	52.067
Conelog®	NP/3.3	14.074		CONELOG NP NE 15.074	52.074	53.100
				CONELOG NP E 16.074	52.074	52.200
	RP/3.8	14.075		CONELOG RP NE 15.075	52.075	53.100
				CONELOG RP E 16.075	52.075	52.200
	WP/5.0	14.076		CONELOG WP NE 15.076	52.076	53.100
				CONELOG WP E 16.076	52.076	52.200

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

 Scan Abutment
on Ti base
53.100 53.200

	Platform	Analogue	Type	Kit	Scan Abutments	
Biohorizns® Internal	NP/3.5	14.087		BH_INTERNAL_NP NE 15.087	52.087	
				BH_INTERNAL_NP E 16.087	52.087	53.100
	RP/4.5	14.088		BH_INTERNAL_RP NE 15.088	52.088	52.200
				BH_INTERNAL_RP E 16.088	52.088	53.100
	WP/5.7	14.089		BH_INTERNAL_WP NE 15.089	52.088	52.200
				BH_INTERNAL_WP E 16.089	52.089	53.100
Biohorizns® External	NP/3.7	14.020		BH_INTERNAL_NP NE 15.051	52.089	52.200
				BH_INTERNAL_NP E 16.051	52.089	53.100
	RP/4.2	14.021		BH_EXTERNAL_RP NE 15.052	50.020	52.020
				BH_EXTERNAL_RP E 16.052	50.020	52.020
	WP/5.2	14.022		BH_EXTERNAL_WP NE 15.053	50.021	52.021
				BH_EXTERNAL_WP E 16.053	50.021	52.021
MIS® Seven	NP/3.7	14.056		MIS_SEV_NP_NE 15.056	50.022	52.022
				MIS_SEV_NP_E 16.056	50.022	52.022
	SP/4.2	14.095		MIS_SEV_SP_NE 15.095 15.295 15.395	50.056	52.056
				MIS_SEV_SP_E 16.095 16.295 16.395	50.056	52.056
	WP/5.2	14.096		MIS_SEV_WP_NE 15.096 15.296 15.396	50.017	52.017
				MIS_SEV_WP_E 16.096 16.296 16.396	50.017	52.017

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

 Scan Abutment
on Ti base
53.100 53.200

	Platform	Analogue	Type	Kit	Scan Abutments
MIS® C1 internal	NP	14.081		MIS C1 NP NE 15.081	52.081 53.100
				MIS C1 NP E 16.081	52.081 52.200
	SP	14.082		MIS C1 SP NE 15.082	52.082 53.100
				MIS C1 SP E 16.082	52.082 52.200
	WP	14.083		MIS C1 WP NE 15.083	52.083 53.100
				MIS C1 WP E 16.083	52.083 52.200
	Multiunit 4.8	14.085		MIS C1 MULTIUNIT WE 15.083	52.085
MIS® V3	NP	14.084		MIS C1 NP NE 15.084	52.084 53.100
				MIS C1 NP E 16.084	52.084 52.200
	SP	14.082		MIS C1 SP NE 15.082	52.083 53.100
				MIS C1 SP E 16.082	52.082 52.200
	Multiunit 4.8	14.085		MIS C1 MULTIUNIT 4.8 15.085	52.085
Megagen AntRidge®		14.057		MGN ANY NE 15.057	50.057 53.100
				MGN ANY E 16.057	50.057 52.200
Anthogyr AXIOM®		14.077		ANTHOGYR Ti base NE 15.077	50.077 53.100
				ANTHOGYR Ti base E 16.077	50.077 52.200
Medentis ICX®		14.080		MED ICX NE 15.080	50.080 53.100
				MED ICX E 15.080	50.080 52.200

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

 Scan Abutment
on Ti base
53.100 53.200

	Platform	Analogue	Type	Kit	Scan Abutments
GlobalD IN-KONE®		14.086		IN_KONE Ti base NE 15.086	52.086 53.100
				IN_KONE Ti base E 16.086	52.086 52.200
				DENTIUM Ti base NE 15.090	52.090 53.100
				DENTIUM Ti base E 15.090	52.090 52.200
Bego Semados®	3.25/3.75	14.091		BEGO Ti base 3.25 NE 15.091	52.091 53.100
				BEGO Ti base 3.25 E 16.091	52.091 52.200
		14.092		BEGO Ti base 4.1 NE 15.091	52.092 53.100
				BEGO Ti base 4.1 E 16.091	52.092 52.200
	4.5	14.093		BEGO Ti base 4.5 NE 15.091	52.093 53.100
				BEGO Ti base 4.5 E 16.091	52.093 52.200
		14.094		BEGO Ti base 5.5 NE 15.094	52.094 53.100
				BEGO Ti base 5.5 E 16.094	52.094 52.200
Biotech Kontakt®		14.078		BIOTECH Ti base NE 15.078	52.078 53.100
				BIOTECH Ti base E 16.078	52.078 52.200
				BTI INTERNAL Base NE 15.079	52.079 53.100
				BTI INTERNAL Base E 15.079	52.079 52.200
Multi-Unit	RP	14.007		Multi_Unit_RP NE 15.007	50.007 52.007
				Multi_Unit_RP NE 15.007/8	50.007 52.007
	WP	14.047		Multi_Unit_RP E 15.054	50.047 52.047

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

 Scan Abutment
on Ti base
53.100 53.200

**Table 1.2 DESS
scan on**

**Uniabutment & Multi-Unit
& Synocta**

	Platform	Analogue	Type	Category	System	Kit	Scan Abutments
Uniabutment® on Astra Tech™ Osseospeed™	Uniabutment® 20°	14.036	●	DESS UNIABUTMENT DIRECT (Uniabutment Scan Abutment)	DESS AST UNIABUTMENT DIRECT (Uniabutment Scan Abutment)	20	50.036 52.036
	Uniabutment® 45°	14.037	●	DESS UNIABUTMENT DIRECT (Uniabutment Scan Abutment)	DESS AST UNIABUTMENT DIRECT (Uniabutment Scan Abutment)	45	50.037 52.037
Astra Tech™ Implant System™ EV	Uniabutment®	14.058	●	DESS UNIABUTMENT DIRECT (Uniabutment Scan Abutment)	DESS AST UNIABUTMENT DIRECT (Uniabutment Scan Abutment)	EV Uniabutment	50.058 52.058

Multi-Unit	RP	14.007	●	DESS MULTIUNIT DIRECT (Multi-Unit Scan Abutment)	DESS AST UNIABUTMENT DIRECT (Uniabutment Scan Abutment)	Multi-Unit	50.036 52.036
	WP	14.047	●	DESS MULTIUNIT DIRECT (Multi-Unit Scan Abutment)	DESS AST UNIABUTMENT DIRECT (Uniabutment Scan Abutment)	Multi-Unit WP	50.037 52.037

Straumann® Soft Tissue level witch synOcta® Abutment	RN on synOcta®	14.010	●	DESS STR SYN DIRECT (Bars and bridges)	DESS STR SYN DIRECT (Bars and bridges)	STR SYN RN 51.010	50.010 52.010
	WN on synOcta®	14.046	●	DESS STR SYN DIRECT (Bars and bridges)	DESS STR SYN DIRECT (Bars and bridges)	STR SYN WN 51.048	50.046 52.037

● non-engaging

⬢ non-engaging

Scan Abutment
for desktop scanning
50.0XX

Scan Abutment
for intra-oral scanning
52.0XX

Scan Abutment
on Ti base
53.100 53.200

Table 1.3 DESS ANGLEBase®

	Platform	Analogue	Type	Kit	Scan Abutments	
Nobel Active®	NP	14.041		NOB_ACT_NP NE 36.041	50.041	52.041
				NOB_ACT_NP E 37.041	50.041	52.041
	RP	14.042		NOB_ACT_RP NE 36.042	50.042	52.042
				NOB_ACT_RP E 37.042	50.042	52.042
Nobel Replace Select™	NP	14.004		NOB_RPL_NP NE 36.004	50.004	52.004
				NOB_RPL_NP E 37.004	50.004	52.004
	RP	14.005		NOB_RPL_RP NE 36.005	50.005	52.005
				NOB_RPL_RP E 37.005	50.005	52.005
	WP	14.006		NOB_RPL_WP NE 36.006	50.006	52.006
				NOB_RPL_WP E 37.006	50.006	52.006
Nobel Brånemark®	NP	14.001		NOB_BRK_NP NE 36.001	50.001	52.001
				NOB_BRK_NP E 37.001	50.001	52.001
	RP	14.002		NOB_BRK_RP NE 36.002	50.002	52.002
				NOB_BRK_RP E 37.002	50.002	52.002
	WP	14.003		NOB_BRK_WP NE 36.003	50.003	52.003
				NOB_BRK_WP E 37.003	50.003	52.003
3i Osseotite®	NP/3.4	14.011		3I_OSS_NP NE 36.011	50.011	52.011
				3I_OSS_NP E 37.011	50.011	52.011
	RP/4.1	14.012		3I_OSS_RP NE 36.012	50.012	52.012
				3I_OSS_RP E 37.012	50.012	52.012
	WP/5.0	14.013		3I_OSS_WP NE 36.013	50.013	52.013
				3I_OSS_WP E 37.013	50.013	52.013

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

	Platform	Analogue	Type	Kit	Scan Abutments	
3i Certain®	NP/3.4	14.014		3I_CER_NP NE 36.014	50.014	52.014 53.100
				3I_CER_NP E 37.014	50.014	52.014 52.200
	RP/4.1	14.015		3I_CER_RP NE 36.015	50.015	52.015 53.100
				3I_CER_RP E 37.015	50.015	52.015 52.200
	WP/5.0	14.016		3I_CER_WP NE 36.016	50.016	52.016 53.100
				3I_CER_NP E 37.016	50.016	52.016 52.200
Straumann® Soft Tissue Level	RN/4.8	14.009		STR_STLV_RN NE 36.046	50.009	52.009
				STR_STLV_RN E 37.009	50.009	52.009
	WN/6.5	14.045		STR_STLV_WP NE 36.045	50.045	52.045
				STR_STLV_WP E 37.045	50.045	52.045
Straumann® Bone Level	NC/3.3	14.043		STR_BLV_NC NE 36.043	50.043	52.043 53.100
				TR_BLV_NC E 37.043	50.043	52.043 52.200
	RC/4.1	14.044		STR_BLV_RC NE 36.044	50.044	52.044 53.100
				TR_BLV_RC E 37.044	50.044	52.044 52.200
Zimmer Screw-Vent®	NP/3.5	14.017		ZIM_SCV_NP NE 36.017	50.017	52.017 53.100
				ZIM_SCV_NP E 37.017	50.017	52.017 52.200
	RP/4.5	14.018		ZIM_SCV_RP NE 36.018	50.018	52.018 53.100
				ZIM_SCV_RP E 37.018	50.018	52.018 52.200
	WP/5.7	14.019		ZIM_SCV_WP NE 36.019	50.019	52.019 53.100
				ZIM_SCV_WP E 37.019	50.019	52.019 52.200

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

 Scan Abutment
on Ti base
53.100 53.200

	Platform	Analogue	Type	Kit	Scan Abutments	
Astra tech implant system™ EV	EV/3.6	14.059		AST_EVO 3.6_NE 36.059	50.059	52.059
				AST_EVO 3.6_E 37.059	53.100	52.200
	EV/4.2	14.060		AST_EVO 4.2_NE 36.060	50.060	52.060
				AST_EVO 4.2_E 37.060	53.100	52.200
	EV/4.8	14.061		AST_EVO 4.8_NE 36.061	50.061	52.061
				AST_EVO 4.8_E 37.061	53.100	52.200
Astra tech osseospeed™	RP	14.024		AST_OSS_RP NE 36.024	50.024	52.024
				AST_OSS_RP E 37.024	53.100	52.200
	WP	14.025		AST_OSS_WP NE 36.025	50.025	52.025
				AST_OSS_WP E 37.025	53.100	52.200
Friadent XIVE®	NP	14.038		DEN_XIV_NP NE 36.038	50.038	52.038
				DEN_XIV_NP E 37.038	53.100	52.200
	RP	14.039		DEN_XIV_RP NE 36.039	50.039	52.039
				DEN_XIV_RP E 37.039	53.100	52.200
	WP	14.040		DEN_XIV_WP NE 36.040	50.040	52.040
				DEN_XIV_WP E 37.040	53.100	52.200

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

 Scan Abutment
on Ti base
53.100 53.200

	Platform	Analogue	Type	Kit	Scan Abutments	
Osstem® TS	MINI	14.072		OSS_TS MINI_NE 36.072	50.072	52.072
				OSS_TS MINI_E 37.072	50.072	52.072
	STD	14.073		OSS_TS MINI_NE 36.073	50.073	52.073
				OSS_TS MINI_E 37.073	50.073	52.073
Camlog®	3.8	14.065		CAM_3.8 NE 36.065	50.065	52.065
				CAM_3.8 E 36.065	50.065	52.065
	4.3	14.066		CAM_4.3 NE 37.066	50.066	52.066
				CAM_4.3 E 37.066	50.066	52.066
Mis® Seven	SP	14.095		MIS_SEV_SP NE 36.017	50.017	52.017
				MIS_SEV_SP E 37.017	50.017	52.017
	WP	14.096		STR_BLV_RC WP NE 36.018	50.018	52.018
				STR_BLV_RC WP E 37.018	50.018	52.018
Megagen AnyRidge®		14.057		MGN ANY NE 36.057	50.017	
				MGN ANY E 37.057	50.017	
Medentis ICX®		14.080		MED ICX NE 36.080 36.180	50.080	
				MED ICX E 37.080 37.180	50.080	
Multi-Unit	RP	14.007		Multi_Unit_RP NE 36.007	50.007	52.007
	WP	14.047		Multi_Unit_RP E 37.054	50.047	52.047

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

 Scan Abutment
on Ti base
53.100 53.200

**Table 1.4 DESS
Pre-milled
Abutment Ti**

	Platform	Analogue	Kit	Scan Abutments	
Nobel Active®	3.0	14.068	NOB_ACT_3.0 NE 61.068	52.068	53.100
	NP	14.041	NOB_ACT_NP NE 61.041 _{ø10mm} 62.041 _{ø14mm}	50.041	52.041
	RP	14.042	NOB_ACT_RP NE 61.042 _{ø10mm} 62.042 _{ø14mm}	50.042	52.042
	WP	14.069	NOB_ACT_WP NE 61.069	53.100	
Nobel Replace Select™	NP	14.004	NOB_RPL_NP NE 61.004	50.004	52.004
					53.100
	RP	14.005	NOB_RPL_RP NE 61.005	50.005	52.005
	WP	14.006	NOB_RPL_WP NE 61.006	53.100	
	6.0	14.097	NOB_RPL_NP 6.0 61.097	52.097	53.100
Nobel Brånemark®	NP	14.001	NOB_BRK_NP NE 61.001	50.001	52.001
	RP	14.002	NOB_BRK_RP NE 61.002	53.100	
	WP	14.003	NOB_BRK_WP NE 61.003	50.003	52.003
3i Osseotite®	NP	14.011	3I_OSS_NP NE 61.011	50.011	52.011
	RP	14.012	3I_OSS_RP NE 61.012	53.100	
	WP	14.013	3I_OSS_WP NE 61.013	50.013	52.013
3i Certain®	NP	14.014	3I_CER_NP NE 61.014	50.014	52.014
	RP	14.015	3I_CER_RP NE 61.015	53.100	
	WP	14.016	3I_CER_WP NE 61.016	50.016	52.016
Straumann® Soft Tissue Level	NNC	14.070	NOB_ACT_WP NE 61.070	52.069	53.100
	RN	14.009	STR_STLV RN NE 61.009 _{ø10mm} 62.009 _{ø14mm}	50.009	52.009
	WN	14.045	STR_STLV WN NE 61.045 _{ø10mm} 62.045 _{ø14mm}	50.045	52.045



non-engaging



non-engaging



Scan Abutment
for desktop scanning
50.OXX



Scan Abutment
for intra-oral scanning
52.OXX



Scan Abutment
on Ti base
53.100

	Platform	Analogue	Kit	Scan Abutments	
Straumann® Bone Level	NC	14.043	STR_BLV_NC NE 61.043 62.043	50.043 53.100	52.043
	RC	14.044	STR_BLV_RC NE 61.044 62.044	50.044 53.100	52.044
Zimmer Screw-Vent®	NP	14.017	ZIM_SCV_NP NE 61.017 62.017	50.017 53.100	52.017
	RP	14.018	ZIM_SCV_RP NE 61.018	50.018 53.100	52.018
	WP	14.019	ZIM_SCV_WP NE 61.019	50.019 53.100	52.019
Astra tech implant system™ EV	EV/3.6	14.059	AST_EVO 3.6_NE 61.059	50.059 53.100	52.059
	EV/4.2	14.060	AST_EVO 4.2_NE 61.060	50.060 53.100	52.060
	EV/4.8	14.061	AST_EVO 4.8_NE 61.061	50.061 53.100	52.061
Astra tech implant Osseospeed™	RP	14.024	AST_OSS_RP NE 61.024	50.024 53.100	52.024
	WP	14.025	AST_OSS_WP NE 61.025	50.025 53.100	52.025
Friadent XIVE®	NP	14.038	DEN_XIV_NP NE 61.038	50.038 53.100	52.038
	RP	14.039	DEN_XIV_RP NE 61.039	50.039 53.100	52.039
	WP	14.040	DEN_XIV_WP NE 61.040	50.040 53.100	52.040
Ankylos® C/X	All platforms, height GH 2,5	14.048 14.049	DEN ANK NE GH2 61.250	50.049 53.100	52.049
Osstem® TS	Mini Yellow	14.072	OSS_TS MINI_NE 61.072	50.072 53.100	52.072
	STD Regular	14.073	OSS_TS STD_NE 61.073	50.073 53.100	52.073



non-engaging



non-engaging



Scan Abutment
for desktop scanning
50.OXX



Scan Abutment
for intra-oral scanning
52.OXX



Scan Abutment
on Ti base
53.100

	<i>Platform</i>	<i>Analogue</i>	<i>Kit</i>	<i>Scan Abutments</i>	
Camlog®	3.3	14.064	CAM_3.3 NE 61.064	52.064 53.100	
	3.8	14.065	CAM_3.8 NE 61.065	50.065 53.100	52.065
	4.3	14.066	CAM_4.3 NE 61.066	50.066 53.100	52.066
	5.0	14.067	CAM_5.0 NE 61.067	52.067 53.100	
Conelog®	NP	14.074	CONOLOG NP NE 61.074	52.074 53.100	
	RP	14.075	CONOLOG RP NE 61.075	52.075 53.100	
	WP	14.076	CONOLOG WP NE 61.076	52.076 53.100	
Biohorizns® External	RP/4.2	14.021	BH_EXTERNAL_RP NE 61.021	50.021 53.100	52.021
	WP/5.2	14.022	BH_EXTERNAL_WP NE 61.022	50.022 52.022	
MIS® Seven	SP	14.095	MIS_SEV_SP_NE 61.017	50.017 53.100	52.017
	WP	14.096	MIS_SEV_WP_NE 61.018	50.018 53.100	52.018
Megagen AntRidge®		14.057	MGN ANY NE 61.057	50.057 53.100	
Dentium Superline™ & Implantium	RP	14.090	DENTIUM Ti base NE 61.090	52.090 53.100	

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

 Scan Abutment
on Ti base
53.100

**Table 1.5 DESS
Pre-milled
Abutment CoCr**

	Platform	Analogue	Kit	Scan Abutments	
Nobel Active®	3.0	14.068	NOB_ACT_3.0 NE 71.068	52.068	53.100
	NP	14.041	NOB_ACT_NP NE 71.041	50.041	52.041
	RP	14.042	NOB_ACT_RP NE 71.042	53.100	
				50.042	52.042
Nobel Replace Select™	WP	14.069	NOB_ACT_WP NE 71.069	53.100	
				52.069	
	NP	14.004	NOB_RPL_NP NE 71.004	50.004	52.004
				53.100	
	RP	14.005	NOB_RPL_RP NE 71.005	50.005	52.005
Nobel Brånemark®	WP	14.006	NOB_RPL_WP NE 71.006	53.100	
				50.006	52.006
	6.0	14.097	NOB_RPL_NP 6.0 71.097	52.097	
				53.100	
	NP	14.001	NOB_BRK_NP NE 71.001	50.001	52.001
3i Osseotite®	RP	14.002	NOB_BRK_RP NE 71.002	53.100	
				50.002	52.002
	WP	14.003	NOB_BRK_WP NE 71.003	50.003	
				52.003	
	NP	14.011	3I_OSS_NP NE 71.011	50.011	52.011
3i Certain®	RP	14.012	3I_OSS_RP NE 71.012	53.100	
				50.012	52.012
	WP	14.013	3I_OSS_WP NE 71.013	50.013	52.013
				53.100	
	NP	14.014	3I_CER_NP NE 71.014	50.014	52.014
Straumann® Soft Tissue Level	RP	14.015	3I_CER_RP NE 71.015	53.100	
				50.015	52.015
	WP	14.016	3I_CER_WP NE 71.016	50.016	52.016
				53.100	
	RN	14.009	STR_STLV RN NE 71.009	50.009	52.009
Straumann® Soft Tissue Level	WN	14.045	STR_STLV WN NE 71.045	50.045	
				52.045	



non-engaging



non-engaging



Scan Abutment
for desktop scanning
50.0XX



Scan Abutment
for intra-oral scanning
52.0XX



Scan Abutment
on Ti base
53.100

	<i>Platform</i>	<i>Analogue</i>	<i>Kit</i>	<i>Scan Abutments</i>	
Straumann® Bone Level	NC	14.043	STR_BLV_NC NE 71.043	50.043 53.100	52.043
	RC	14.044	STR_BLV_RC NE 71.044	50.044 53.100	52.044
Zimmer Screw-Vent®	NP	14.017	ZIM_SCV_NP NE 71.017	50.017 53.100	52.017
	RP	14.018	ZIM_SCV_RP NE 71.018	50.018 53.100	52.018
	WP	14.019	ZIM_SCV_WP NE 71.019	50.019 53.100	52.019
Astra tech implant system™ EV	EV/3.6	14.059	AST_EVO 3.6_NE 71.059	50.059 53.100	52.059
	EV/4.2	14.060	AST_EVO 4.2_NE 71.060	50.060 53.100	52.060
	EV/4.8	14.061	AST_EVO 4.8_NE 71.061	50.061 53.100	52.061
Astra tech implant Osseospeed™	RP	14.024	AST_OSS_RP NE 71.024	50.024 53.100	52.024
	WP	14.025	AST_OSS_WP NE 71.025	50.025 53.100	52.025
Friadent XIVE®	NP	14.038	DEN_XIV_NP NE 71.038	50.038 53.100	52.038
	RP	14.039	DEN_XIV_RP NE 71.039	50.039 53.100	52.039
	WP	14.040	DEN_XIV_WP NE 71.040	50.040 53.100	52.040
Ankylos® C/X	All platforms, height GH 2,5	14.048 14.049	DEN ANK NE GH2 71.250	50.049 53.100	52.049
Osstem® TS	Mini Yellow	14.072	OSS_TS MINI_NE 71.072	50.072 53.100	52.072
	STD Regular	14.073	OSS_TS STD_NE 71.072	50.073 53.100	52.073



non-engaging



non-engaging



Scan Abutment
for desktop scanning
50.0XX



Scan Abutment
for intra-oral scanning
52.0XX



Scan Abutment
on Ti base
53.100

	<i>Platform</i>	<i>Analogue</i>	<i>Kit</i>	<i>Scan Abutments</i>	
Camlog®	3.3	14.064	CAM_3.3 NE 71.064	52.064 53.100	
	3.8	14.065	CAM_3.8 NE 71.065	50.065 53.100	52.065
	4.3	14.066	CAM_4.3 NE 71.066	50.066 53.100	52.066
	5.0	14.067	CAM_5.0 NE 71.067	52.067 53.100	
Conelog®	NP	14.074	CONELOG NP NE 71.074	52.074 53.100	
	RP	14.075	CONELOG RP NE 71.075	52.075 53.100	
	WP	14.076	CONELOG WP NE 71.076	52.076 53.100	
Biohorizns® External	RP/4.2	14.021	BH_EXTERNAL_RP NE 671.021	50.021 53.100	52.021
	WP/5.2	14.022	BH_EXTERNAL_WP NE 71.022	50.022 52.022	
MIS® Seven	SP	14.095	MIS_SEV_SP_NE 71.017	50.017 53.100	52.017
	WP	14.096	MIS_SEV_WP_NE 71.018	50.018 53.100	52.018
Megagen AntRidge®		14.057	MGN ANY NE 71.057	50.057 53.100	
Dentium Superline™ & Implantium	RP	14.090	DENTIUM Ti base NE 71.090	52.090 53.100	

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
50.0XX

 Scan Abutment
for intra-oral scanning
52.0XX

 Scan Abutment
on Ti base
53.100

All DESS products are manufactured according to ISO 9001:2015 and ISO 13485:2016 with CE marking for all products Class IIb

All DESS® AURUMBas® screws and screwdrivers are based on the standard ISO 10664

3i® and 3i Certain®, Miniplant®, MicroMiniplant™ and Osseotite® are registered trademarks of 3i Implant Innovations, USA

ASTRA TECH™, OsseoSpeed™, ANKYLOS® C/X, XiVE®, FRIADENT®, ASTRATECH implant system™ EV and UniAbutment

are registered trademarks of DENTSPLY SIRONA Implants

BioHorizons® is a registered trademark of BioHorizons, INC.

Brånemark System®, NobelReplace® and Replace Select™, NobelActive™, NobelReplace® CC, NobelSpeedy™ and Multi-Unit are

registered trademarks of Nobel Biocare AB, Sweden

ITF®, Straumann® and synOcta® are registered trademarks of Straumann Holding AG, Switzerland

MIS® is a registered trademark of MIS Implant Technologies Ltd.

Zimmer Screw-Verband SwissPlus® are registered trademarks of Zimmer Dental INC, USA

Dyna® is a registered trademark of Dyna Dental Engineering BV

Locator® is a registered trademark of ZEST IP Holdings

Camlog® and Conelog® are registered trademarks of CAMLOG Biotechnologies AG

Megagen® is a registered trademark of MEGAGEN IMPLANT CO. LTD

TORX® is a trademark of ACUMENT INTELLECTUAL PROPERTIES, LLC

Bego Semado® is a registered trademark of BEGO GmbH & Co. KG

Anthogyr® is a registered trademark of Anthogyr SAS

Medentis ICX is a registered trademark of medentis medical GmbH

Osstem® is a registered trademark of OSSTEM IMPLANT CO. LTD

Dentium® is a trademark of Dentium Co. Ltd.

GOBALD® is a trademark of Global D SASU.

Biotech Kontakt is a trademark Biotech Dental SAS.

BTI® is a trademark BTI Biotechnology Institute.

DESS products can be used only by licensed dentists, physicians and labs.

Please refer ever to the instructions of use available on the website

www.dess-usa.com

or supplied with the product.

Some products may not be available in your country. Some of the products are cleared for US and Canada.

Please contact your local distributor for more details.

